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
Thank you for purchasing a **pH combination electrode**.

For the best results with your electrode:

- **read** this user manual carefully,
- **comply** with the precautions for use.

This electrode is used to measure the pH of a solution.

 The CE marking indicates compliance with the European directives, particularly LV and EMC.

 The crossed-out dustbin means that, in the European Union, the product is the subject of selective disposal as per directive DEEE 2002/96/EC: this equipment must not be treated as household waste.

USE

- Remove the storage reservoir cap.
- If you notice the presence of an air bubble at the level of the glass membrane or diaphragm, remove it by shaking the electrode.
- Connect the electrode to the measuring instrument.
- Open the opening tab.
- Calibrate the electrode with buffer solutions, rinsing with demineralized water after each measurement.
- Immerse the tip of the electrode in the solution to be measured, while keeping it vertical and making sure that the membrane and diaphragm are immersed, but without exceeding the level of the internal electrolyte.
- Wait for the measurement to stabilize.
- Rinse the electrode with demineralized water after each measurement.
- Close the tab again when you have finished measuring.


STATE AT DELIVERY

The pH electrode is delivered in a cardboard box with:

- a shim,
- a storage reservoir cap,
- a bilingual user's manual,
- a multilingual accompanying datasheet.

For accessories and replacement parts, please visit our website:

www.chauvin-arnoux.com

The presence of white crystalline deposits on the electrode is a normal phenomenon due to evaporation of the storage solution in the reservoir cap or the internal electrolyte. To remove this deposit, simply rinse the electrode with demineralized water and then wipe it off with a clean, lint-free cloth without touching the glass membrane.

Do not use aggressive and/or abrasive products (hydrofluoric acid, acetone, etc.). The glass membrane must not be scratched during cleaning.

Contamination of the diaphragm is the most frequent cause of measurement problems.

- Remove any impurities deposited on the glass surface of the electrode and on the outside of the diaphragm with demineralized water or with a soft cloth if the deposit is significant.
- If rinsing is not sufficient, it is possible to use various chemical cleaning solutions, depending on the type of contamination :

Type of clogging	Electrodes concerned	Cleaning solution
By silver sulphide	Electrode filled with a KCl electrolyte	Solution of thiourea in hydrochloric acid
By silver chloride		Concentrated ammoniac solution
By proteins	Any electrodes with an electrolytic compartment	Solution of pepsin + HCl
Other types of clogging		Ultrasonic bath with water or HCl 0.1 mol/L

MAINTENANCE

Keep your electrode in a vertical position in a reservoir cap containing an electrolytic solution identical to the electrode's internal electrolyte (KCl 1 mol/L in most of our electrodes). The level of the storage electrolyte must be lower than that of the internal electrolyte with which the electrode is filled.

Electrode reference	Filling electrolyte	Storage solution
BRV1A/BRV1H BRV4A/BRV4H BRV5A/BRV5H BRV15H BRV22A/BRV22H	KCl 1 mol/L	KCl 1 mol/L storage solution
BRSb1		
DRV2A/DRV2H		
LRV22H		
XRV1A/XRV1H XRV22H XRVST1H	Polymer	Glycerated KCl 1 mol/L solution
LRV6H		
BRV45H	Acetic acid	KCl 1 mol/L storage solution

For electrodes with liquid electrolytes, close the tab of the filling hole before storing.

Do not store an electrode in dry conditions or in demineralized water. The membrane may be affected and the life span of the electrode may be reduced.

Regularly check the condition of your electrode, the cleanness of the diaphragm and glass membrane and the level of the electrolyte.

If the level of the electrolyte in the electrode seems insufficient or contaminated, you will have to top it up or empty it and refill it with a new electrolyte solution suitable for your electrode.

LIFE SPAN AND WARRANTY

pH electrodes are wear parts. The life span of your electrode depends on the conditions in which it is used and the maintenance which you carry out on it.

The warranty does not apply if the glass is broken.

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