



Operation Manual

Cubic Vacuum Desiccators

**Models: VDC-11, 21, 31, 41,
11U, 21U, 31U, 41U**

Manual No. D4001L002, Version 0.0

Table of Contents

1	Safety Advice	1
1.1	Operation Manual	1
1.2	DANGER/CAUTION/NOTICE Alerts	1
1.3	Danger Alerts.....	1
1.4	Caution Alerts.....	2
2	Product Description.....	3
2.1	Introduction.....	3
2.2	Characteristics	4
2.3	Layout.....	6
3	Unpacking and Installation.....	8
3.1	Checking Shipment Damage and Unpacking.....	8
3.2	Checking the Contents of the Package.....	8
3.3	Installation.....	8
3.3.1	Installation Environment	8
3.3.2	Installation location	9
4	Operation.....	10
4.1	Precautions	10
4.2	Operation of the 3-way Valve	10
4.3	Vacuum Formation.....	11
4.3	Multi-stacking of the Units.....	12
4.4	Storage.....	12
5	Maintenance	13
6	Troubleshooting.....	14
7	Accessories	15
8	Appendix.....	16
8.1	Technical Specifications	16
8.2	Disposal of the Unit.....	17
8.3	Warranty and Disclaimer	17
8.4	After-sales Service and Customer Assistance.....	18

1 Safety Advice




1.1 Operation Manual

This manual contains important safety and operation information. You must carefully read, understand, and follow all the instructions in this manual prior to operating this instrument. Keep this manual in a safe place nearby for reference and make it easily available to all users.

1.2 DANGER/CAUTION/NOTICE Alerts


This manual highlights DANGER/CAUTION/NOTICE alerts to prevent injury or property damage and also to achieve optimum performance of your instrument.

These alerts are classified into three types in this manual depending on the importance and the risk levels as described below:

Symbols	Meaning
	Ignoring this alert could cause serious or even fatal injury.
	Ignoring this alert could cause serious injury or property damage.
	Ignoring this alert could cause operational problems.

Ensure to read the danger and caution alerts in 1.3 and 1.4 carefully before use.

1.3 Danger Alerts


Ignoring the following warnings could cause serious or even fatal injury..

- ◆ **Never install or use this unit in explosive atmospheres.**
- ◆ **Never install this unit near to hazardous or flammable substances.**
- ◆ **Never disassemble, repair, or modify this unit on your own. Doing so will void your warranty and may result in injuries or product damages.**

1.4 Caution Alerts

 CAUTION
Ignoring the following cautions could cause serious injury or property damage.

- ◆ **Carefully read all the warning labels before using this unit. And do not remove or damage the warning labels.**
- ◆ **Do not move this unit during operation. When moving it, use the molded-in handles.**
- ◆ **Do not use this unit in wet or moist atmosphere or where water leakage is expected.**
- ◆ **Do not use this unit in contaminated atmosphere or where metallic dust exists.**
- ◆ **Do not expose this unit to any heat sources including direct sunlight.**
- ◆ **Do not allow moisture, organic solvent, dust, or corrosive gases to get into this unit.**
- ◆ **If you observe a strange smoke, odor or noise from this unit, stop using it. After the smoke or odor disappears, contact your dealer or Jeio Tech if any repair is required.**
- ◆ **Do not use chlorine bleach, ammonia-based cleaners, abrasives, ammonia, or metal scouring pads. Wipe with a soft dampened cloth or a sponge soaked in water or diluted neutral detergent.**

2 Product Description

2.1 Introduction

Congratulations of your purchase of our high-quality cubic vacuum desiccator which is specially designed to provide maximum benefit for your investment with respect to performance, safety, ease of use, and durability.

This unit is a sturdily constructed 360° transparent vacuum desiccator ideal for storing, drying, cooling, and degassing various materials. Vacuum environment created by this unit providing a very low humidity, contamination-free, corrosion-free, and dust-free environment is ideal for:

- storing moisture-sensitive materials such as anhydrides and other hygroscopic compounds as well as materials which are highly reactive with atmosphere,
- removing residual solvents or air bubbles,
- cultivating anaerobic microbes,
- minimizing the reabsorption of humidity and the exposure to air during the drying or the cooling stages of biological/medical/food/materials experiments,
- conducting experiments at low pressure,
- drying at low pressure to minimize the heat effects.

Currently eight models of cubic vacuum desiccators are available from Jeio Tech (VDC-11/11U, 21/21U, 31/31U, and 41/41U) with the internal volume ranging from 11 liter to 41 liter. Note that UV blocking models (VDC-11U, 21U, 31U, and 41U) are offered to minimize damages or discoloration of light-sensitive samples by using UV-blocking resin which completely blocks all ranges of UV and the blue spectrum of the visible light.

2.2 Characteristics

◆ Key Vacuuming Features

(1) Outstanding Vacuum Sustainment

This unit maintains a 1 Torr (133 Pa) vacuum for over 72 hours at room temperature and stays gastight without a vacuum allowing reliable experiments over extended periods.

(2) Integral Molding of the Polycarbonate Body

Unlike conventional vacuum desiccators, the polycarbonate cubic body (except the door part) of this unit is integrally molded as a single piece for maximum sealing performance. This unique design without requiring any adhesive for body construction allows virtually no leakage.

Shock-resistant, high-tensile strength, and transparent polycarbonate body provides extreme durability as well as 360° visual observation from outside.

(3) Viton® 3-way Valve

This unit comes standard with a Viton® 3-way valve offering consistent and uniform vacuum draw, vacuum release, or gas exchange without the inconvenience of connecting and disconnecting pump hoses to the unit.

If needed, up to 4 valves can be installed using the built-in ports at the base of the unit and the positioning of the valve(s) can be also easily changed.

(4) Airtight Door Sealing

The door is located in the front for convenient access. Equipped with a special door latch which is easily closed at normal pressure and a greaseless vacuum seal using high-quality silicon gasket, this unit enables airtight sealing between the door and the body and quick vacuum formation.

◆ **Key Design Features**

(1) Multi-stackable Design

The VDC models come in a variety of capacities to suit the user's needs with the internal volume ranging from 11 liter to 41 liter.

With the same foot print, different models can be stacked together in any convenient fashion allowing for lab space saving.

(2) Efficient Use of the Internal Space

Compared to round desiccators, cubic desiccators provide larger usable volumes with the same footprint. In addition, the internal space can be more efficiently utilized by using the standard perforated shelves designed for efficient airflow.

Two or three perforated shelves are provided and, if needed, up to four or six shelves can be used for efficient storage of the samples. (Refer to Product Checklist or Specifications for details.)


(3) Easy Vacuum Level Checking

With an built-in vacuum gauge installed on the front door, the vacuum level checking can be done conveniently.

(4) Desiccant Tray

This unit comes standard with a desiccant tray.

2.3 Layout



- 1. Vacuum Gauge**
- 2. Door**
- 3. Door Latch**
- 4. 3-way Valve**
- 5. Non Slip Feet**
- 6. Groove Handle**
- 7. Valve Port Plug**
- 8. Door Seal**
- 9. Shelf Rail**
- 10. Perforated Shelves**
- 11. Desiccant Tray**
- 12. Stacking Groove**
- 13. Air Inlet/Outlet Nozzle**









3 Unpacking and Installation

3.1 Checking Shipment Damage and Unpacking

Upon receiving the instrument, check to ensure that no damage has occurred during shipment. It is important that any damage that occurred during shipment is detected before unpacking. If you find such damage, the carrier must be notified immediately.

3.2 Checking the Contents of the Package

After unpacking, check to ensure that all the parts and accessories described below are included in the package. If not, contact your dealer immediately.

Item		Figure	Quantity
Main Body			1
Vacuum Gauge			1
3-way Valve			1
Port Plug	VCD-11/11U		1
	VDC-21/21U, 31/31U, 1/41U		3
Door Seal			1
Desiccant Tray			1
Perforated Shelves	VCD-11/11U, 21/21U		2
	VDC-31/31U, 41/41U		3
Operation Manual			1

3.3 Installation

3.3.1 Installation Environment

NOTICE

Avoid exposure to any heat sources including direct sunlight.

This unit is designed for indoor use in laboratory environments. The conditions for proper installation are as follows:

Permissible Ambient temperature	5 to 40°C (41 to 104°F)
Permissible Relative humidity	below 80%
Permissible Altitude	0 to 2,000 m (6,562 ft)

3.3.2 Installation location

NOTICE




Install the unit on a level surface.

4 Operation

4.1 Precautions

⚠ CAUTION
<ul style="list-style-type: none"> ● Do not form a vacuum when airtight containers are inside the unit. ● Note that this unit cannot be autoclaved. ● When attaching other equipments or additional 3-way valves after removing the port plugs, use proper O-rings and nuts for tight sealing. ● When putting high-temperature samples inside the unit, cool them thoroughly in advance or make sure that they are not in direct contact with the interior of the unit to avoid interior damage.

4.2 Operation of the 3-way Valve

	Front Nozzle	Rear Nozzle
	Open	Closed
	Closed	Open
	Closed	Closed

4.3 Vacuum Formation

The procedures for vacuum formation are as follows:

- (1) When using the unit for drying purposes, fill the tray with desiccants and put it inside prior to operation.
- (2) After putting the samples inside the unit, firmly close the door by locking the two latches properly. You will hear a click sound when the latch is properly engaged.
- (3) Connect the vacuum pump hose to one of the two nozzles of the 3-way valve.
- (4) After activating the vacuum pump, turn the valve handle towards the nozzle to open it and read the vacuum gauge to check whether the vacuum is forming properly.
- (5) When the desired vacuum level is achieved, turn the valve handle until it is perpendicular to the nozzle to close it and check the vacuum gauge whether there is any leakage.
- (6) When releasing the vacuum, slowly turn the valve handle towards the opposite nozzle to return to the normal pressure.

CAUTION


- Do not impact or drop the unit when a vacuum is formed inside.
- When forming or releasing a vacuum, turn the valve handle slowly to avoid rapid pressure change. If not, the stored samples can be damaged or dispersed.
- When forming a vacuum, install a cold trap in front of the vacuum pump. Moisture and chemicals from the samples may result in damages to the pump.
- Do not let the body of the unit into direct contact with acetone, benzene, toluene, chloroform, cresol, sodium hydroxide, highly concentrated nitric or sulfuric acid, acetic acids, or strong chlorine-based solvents.
- Ensure that the door seal is not damaged by strong acid or base. Damaged seal can cause malfunction of the unit.
- When a vacuum is formed inside the unit, do not apply excessive force to the 3-way valves, the vacuum gauge, or any of their vicinities. Any damages

caused by external forces can cause malfunction of the unit.

- When the unit is not forming a vacuum, keep the valve handle perpendicular to the nozzles.
- When relocating the unit, use the built-in groove handles instead of holding on to other protruded parts such as the vacuum gauge or 3-way valves.

4.3 Multi-stacking of the Units

All VDC models are designed to be easily multi-stacked together by locking the feet of the unit into the stacking grooves located on the top surface of another unit as shown below:



4.4 Storage

If this unit is not used for an extended period of time, clean the instrument with soft cloth and store it in dry place after packing it properly.

5 Maintenance

 **CAUTION**

- If the unit is contaminated, wear chemicals-resistant gloves before cleaning.
- Regularly check the 3-way valves, the vacuum gauge, or any of their vicinities for any damage.
- Regularly check the door seal as well as the surface where the door and the body are in contact and keep them always clean and undamaged.
- Make sure to clean the door seal regularly using diluted neutral detergent. For thorough cleaning, detach the seal from the door.
- Do not use chlorine bleach, ammonia-based cleaners, abrasives, organic solvents, or metal scouring pads when cleaning. Use a soft cloth all the time.

6 Troubleshooting





Problem	Cause	Corrective Action
Vacuum Malfunction	Contaminated surface between the door and the body	Thoroughly clean the surface with a soft cloth using diluted neutral detergent.
	Damaged door seal	Replace the seal.
	Unsuitable room temperature	Check the temperature whether it is within the operating range.
	Misplaced or damaged 3-way valve	Check whether the valve is properly set. If damaged, replace it with a new one.
	Damaged vacuum gauge	Replace the gauge.

7 Accessories

Component	Cat. No.	Description	
Vacuum Gauge	EDA4214	0 to -0.1 (Mpa)	
Door Seal	VDC0032	VDC-11(U)	Silicone, 220 x 250 x 6.5
	VDC0033	VDC-21(U)	Silicone, 310 x 260 x 6.5
	VDC0015	VDC-31(U)	Silicone, 310 x 360 x 6.5
	VDC0034	VDC-41(U)	Silicone, 310 x 460 x 6.5
Desiccant Tray	VDC0037	VDC-11(U)	PP, 158 x 24 x 194
	VDC0005	VDC-21(U)/31(U)/41(U)	PP, 268 x 218 x 26
Perforated Shelves	VDC0035	VDC-11(U)	PC, 202 x 6.5 x 204
	VDC0003	VDC-21(U)/31(U)/41(U)	PC, 289 x 11 x 304
3-way Valve	AAAD4501	PP Ø9.5 x 25	
Silicagel (20g)	VDR0004	For all VDC models	
Silicagel (40g)	VDR0005	For all VDC models	

8 Appendix

8.1 Technical Specifications

Model		VDC-11(U)	VDC-21(U)	VDC-31(U)	VDC-41(U)
Figure					
Internal volume (L)		11	23	35	45
Vacuum	Gauge range (MPa)	0 ~ -0.1			
	Max. Permissible	1.33 x 10 ⁻⁴ (1Torr)			
	Nozzle diameter (mm)	Ø 9.5			
Material	Body, Shelves	Polycarbonate, transparent			
	3-way Valve, Latch, Desiccant tray	Polypropylene			
	Vacuum Seal, Foot	Silicone Rubber			
Dimensions (mm)	Overall (W x D x H)	322x285x271	420x397x281	420x397x381	420x397x491
	Door opening (W x H)	208x238	295x245	295x345	295x445
	Shelf (W x D x H)	202x204x6.5	289x304x11		
	Desiccant tray(W x D xH)	158x194x24	218x268x268		
Net weight (kg)		4.2	8.2	10.8	12.2
Shelves	Number of shelves (standard/maximum)	2/4	2/4	3/5	3/6
	Max. load/shelf (kg)	3	5	5	5
Model		VDC-11U	VDC-21U	VDC-31U	VDC-41U
Color		Amber			
UV Protection (range)		UV A,B,C 100% (200~450nm)			

8.2 Disposal of the Unit

Disposing of the unit must be done in an environmentally responsible way if it has been potentially exposed to bio-agents or radioactive samples. Failure to follow stringent requirements for appropriate disposal may lead to actions against you and your organization.

- (1) First, check with your laboratory or organization to ensure that you are following all the policies and procedures for disposal of laboratory equipments and wastes.
- (2) If not possible, contact your local governing body for regulations regarding disposal of laboratory equipments. Jeio Tech highly recommends you to find a local service provider that can properly dispose of your unit.

8.3 Warranty and Disclaimer

Jeio Tech warrants to the original consumer that the products it manufactures will be free of defects and perform as specified for a period of two year from the purchase date if the unit was used according to the instructions in this operation manual and only with accessories or consumable parts approved by Jeio Tech.

Any of Jeio Tech's products which, under normal operating conditions, proves defective in material or in workmanship, within two year from the purchase date will be repaired or replaced free of charge.

Our warranty excludes accessories or consumable parts and does not apply to faults resulting from abuse, misuse, alteration, insufficient care, maintenance contrary to the instructions in this operation manual.

Warranty claims can be processed faster if the unit's serial number, detailed description of the problem, and its operating conditions are provided. Consumers will be required to submit the original receipt as proof of purchase date to support a warranty claim.

JEIO TECH HEREBY DISCLAIMS ALL OTHER WARRANTIES, WHETHER WRITTEN OR ORAL, EXPRESS OR IMPLIED BY LAW OR OTHERWISE, INCLUDING WITHOUT LIMITATION, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.

CUSTOMER'S SOLE REMEDY FOR ANY DEFECTIVE PRODUCT WILL BE AS STATED ABOVE, AND IN NO EVENT WILL JEIO TECH BE LIABLE FOR INCIDENTAL, CONSEQUENTIAL, SPECIAL OR INDIRECT DAMAGES IN CONNECTION WITH THE PRODUCT.

8.4 After-sales Service and Customer Assistance

Our after-sales service responds to your questions concerning maintenance and repair of your unit as well as spare parts. In case of a warranty claim, repair or purchase of replacement parts or in case of queries or other problems, please contact your local dealer or Jeio Tech representative.

International Sales Head Office (Korea)

#1005, Byucksan Digital Valley 6-Cha, 481-4 Gasan-Dong, Geumcheon-Gu, Seoul 153-704, Korea

Tel: +82 2 2627 3816 E-mail: overseas@jeiotech.com

FAX: +82 2 3143 1824

The Americas (U.S.A. Branch)

1-A Gill St. Woburn, MA 01801, U.S.A.

Tel: +1 781 376 0700 E-mail: info@jeiotech.com

FAX: +1 781 376 0704

Europe (U.K. Branch)

Unit 3, Tower Industrial Park, Chalgrove, Oxfordshire, OX44 7XZ, United Kingdom

Tel: +44 1865 400321 E-mail: labcompanion@medlinescientific.com

FAX: +44 1865 400736

China (Shanghai Branch)

A-2113 Oriental International Plaza, 85 LouShanGuan Rd, Changning District, Shanghai, China 200336

Tel: +86-21-5108-9161 +86 21 3251 1086 E-mail: longjuncao@jeiotech.com

FAX: +86 21 3251 1083

South East Asia (Malaysia Branch)

No 57-59, Jalan Adenium 2G/6, Pusat Perniagaan Adenium, 48300 Bandar Bukit Beruntung, Selangor Darul Ehsan, Malaysia

Tel: +60 3 6028 5833 +60 3 6028 5825 E-mail: labcomp@streamyx.com

FAX: +60 3 6028 5822

◇ The contents of this manual can be changed or upgraded without prior notice.
◇ The copyright of this manual is reserved by Jeio Tech.