

Laboratory conductivity meter

CA 10141E

**Laboratory
conductivity meter**



- ✓ Compatible with **Regressi**, **Graph2D** and **DLT** via **USB port**
- ✓ Compatible with **CAEx** interfaces via the **analogue outputs**
- ✓ Storage of **100,000** measurements
 - ✓ Ideal for **fixed and mobile measurements**
 - ✓ **Rugged and modern**

TARGET MARKET

Education: priority market

Secondary education: mainly high schools. For the general and technological/scientific sections (e.g. STL, STI2D, ST2S). May also be proposed in middle schools.

Higher education: Universities (Bachelor's, Master), IUT (BUT) technological colleges, high schools (CPGE and BTS), Engineering schools. In general terms, any establishments proposing a scientific syllabus such as Chemistry, Physics, Life Sciences, Biology, Biotechnologies, Engineering, Health, Pharmacy, Environment, Agri-food, etc.

Conductivity meters are used during **Practical Exercises**. In these Practicals, the students perform scientific experiments.

Analogue outputs: enable connection to CA Ex (Computer Assisted Experimentation) interfaces. A CAEx interface is a unit connected to a PC and to sensors (pH-meter, conductivity meter, thermometer, oxygen analyser, etc.). These allow users to acquire and process the experimental data. They are sometimes accompanied by specific software which create a virtual laboratory for the student.

Regressi and Graph2D: educational software. In particular, they can be used to acquire and process experimental data. The pH-meter can be hooked up to these software products (1 at a time!) which record the values measured. Students can then process these values (as a curve, for example).

Data Logger Transfer (DLT): can be used to recover the saved measurements, configure the instrument and program the recordings.

Any activities not requiring the industrial range's IP 67 protection: secondary market

Academic research laboratories
Public and industrial R&D laboratories
Regulatory test laboratories

Sectors: agri-food, environment, agriculture, drinking water, chemicals industry, cosmetics, biotechnology, etc.

TARGET USERS

Teachers, students, researchers, engineers, technicians, chemists, etc.

COMPETITIVE ADVANTAGES

- ✓ Compatible with **Regressi** and **Graph2D** via **USB port**.
- ✓ Compatible with **CAEx** interfaces via the **analogue outputs**.
- ✓ Compatible with **Data Logger Transfer** (programmable recordings).
- ✓ Automatic or manual storage of **100,000 time/date-stamped measurements**.
- ✓ Micro USB port for **transfer onto computer**.
- ✓ Extra-wide multi-display **LCD screen**: clear display for easy reading of values.

- ✓ **Guided calibration** (1 point) with a **customizable** list of buffer solutions.
- ✓ **Shockproof** sheath, **backlighting**.
- ✓ Made in **France**, **eco-design**.

SPECIFICATIONS

Conductivity	Measurement range	0.050 μ S/cm to 200.0 mS/cm
	Resolution (R)	1 nS/cm to 100 μ S/cm (depends on range)
	Intrinsic uncertainty*	$\pm 1\% r^{**} \pm R$
TDS	Measurement range	0.001 mg/l to 200.0 g/l
	Resolution (R)	1 μ g/l to 100 mg/l (depends on range)
	Intrinsic uncertainty	$\pm 1\% r \pm R$
Resistivity	Measurement range	2 Ω .cm to 4999 k Ω .cm
	Resolution (R)	1 m Ω .cm to 1 k Ω .cm (depends on range)
	Intrinsic uncertainty	$\pm 1\% r \pm R$
Salinity	Measurement range	2.0 to 42.0 psu
	Resolution (R)	0.1 psu
	Intrinsic uncertainty	$\pm 0.5\% r \pm R$
Temperature	Measurement range	-10 to 120 °C / 14 to 248.0 °F
	Resolution (R)	0.1 °C / 0.1 °F
	Intrinsic uncertainty	$\pm 0.4^\circ\text{C} / \pm 0.7^\circ\text{F}$
	Tref available	20/25°C (68/77°F)
Calibration	1 point, 6 predefined conductivity references modifiable by the user	
T°C compensation	Automatic (ATC) or manual (MTC), -10 to 120°C / 14 to 248 °F	

Temperature correction	Linear, non-linear, no correction
Probe	NOT SUPPLIED
Data storage	100,000 measurements
Batteries / Battery life	4 x 1.5 V AA or LR6 alkaline batteries / ~200h use
Sensors input	BNC (conductivity cell) Jack (Pt1000 temperature probe)
Communication interface	Type-B micro-USB (also used for power supply) Analogue outputs for conductivity/TDS/salinity/resistivity and temperature (3 x 4 mm banana)
Dimensions (with sheath)	211x127x54 mm
Weight	600 g (with batteries)
Environmental conditions	Storage: -20 to +70°C Operation: -10 to +55 °C
Warranty	3 years

*The intrinsic uncertainties on the measurements are indicated for the instrument alone.

**r : % of reading

STATE AT DELIVERY

The **CA 10141E conductivity meter** is delivered in a cardboard box with:

- 1) 1 protective sheath mounted on the instrument
- 2) 4 x AA or LR6 alkaline batteries
- 3) 1 USB-mains adapter
- 4) 1 USB – micro USB cable
- 5) Bilingual Quick Start Guide
- 6) Verification certificate

The cells and probes are not supplied with the instrument.

TO ORDER

CA 10141E laboratory conductivity meter.....P01710021

Accessories

Glass/platinum conductivity cellBCP4-BNC
PVC/platinum conductivity cell.....XCP4-BNC
Temperature sensor (Pt1000).....P01710070
1408 μ S/cm conductivity reference solutionP01700118*
147 μ S/cm conductivity reference solution, 125 mLP01700117*
12.85 mS/cm conductivity reference solution, 125 mL.....P01700119**
KCl 1mol/L conductivity reference solution, 125 mL.....P01700116
Set of 3 plastic beakers.....P01710056

Replacement parts

Shockproof sheath.....P01710050
USB- μ USB cable & mains adapter ...P01654023

*delivered with a quality certificate guaranteeing traceability to the NIST (National Institute of Standards and Technology, USA) reference materials

** delivered with a quality certificate guaranteeing traceability to the OIML (Organisation Internationale de Métrologie Légale) reference materials