## **WIM**

# **Mehrkammer-Inkubator**

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# **Multi-Room Incubator**

## **Operating Manual**

gilt für / written for  
WIM -4/ -R4/ -RL4

**witeg Labortechnik GmbH**

(English manual at page 11)

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# 1 Introduction

Thank you for purchasing our **WIM Multi-Room Incubator**. You did a good choice buying a **Wisd** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a save place for future reference and hand it out if you consign the device to a third person.

**Witeg** reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

## 2 Safety Instructions

### 2.1 Description of symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

**ATTENTION to high temperatures**

Instructions to avoid injuries caused by high temperatures.

### 2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.

- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- Pay attention if you open the door, hot air is coming out.
- Do not put heat-sensitive items at the device.
- Warn all attendees if you leave the device unattended.
- Do not clean the scorching device.

- Repairs and amendments should only be made by **witeg** authorized people. Contraventions exclude the warranty.
- Make sure, that the safety switch is adjusted to 20-30°C over max. temperature. Otherwise the heating maybe doesn't work.

- Use this device only if you read the whole safety instructions.
- This device should only be used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reach the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- All circuits and electrotechnical components are **witeg** owned technology. In fact of this it is necessary that only **witeg** authorized personal repair the device.
- Pull out the power cable if you didn't want to use the device for a longer time.

- It is important that the device stands up on a horizontal and firm subsurface to avoid vibrations.
- If you set the device in a site with vibrations, make sure that the device can not slip away.
- The device could get damaged if the vibrations are too strong.
- To avoid heat congestion or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.

- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superiorly to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- Do not use this device in area with potentially explosive atmosphere.

### 3 Packing contents

quantum	denomination
1	<b>WIM Multi-Room Incubator</b>
8	<b>PTFE-coated shelve</b>
1	<b>operating manual</b>

### 4 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your customer service.
4. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.
5. Make sure that you follow all safety instructions.
6. Connect the device with the power outlet.
7. Switch the device on.
8. Set your values and start the process (look at job description).

### 5 Product description

#### 5.1 Key features

The **Wisdom WIM Multi-Room Incubator** ensures a constant temperature which makes it suitable for microorganism culture, germination and animal and plant cell culture from 0°C up to 70 °C. The device was developed for simplest usage and a very long life.

Key features are as follows:


- Simplest control by Jog-Shuttle knob.
- Fast commissioning.
- Long life by simply and solid construction.
- Time and temperature showed on illuminated LCD.
- Storage functions for adjusted time and temperature values.
- Turn-off and turn-on delay from one minute up to over 99 hours.
- Signal sound after end of delay time.
- Easy to justify.
- Key lock against inadvertent changing of set temperature and time values.
- Height-adjustable PTFE-coated steel shelves included.


- Over temperature and over current protection.
- Inner room made of stainless steel.
- Glass door for process viewing without interruption.

## 5.2 General survey

The functionality of the described components is the same at each type of the **WIM**-series. Only the position differs from type to type.

## 5.3 Operating surface

## 5.4 Term description

**SV:** set value (upper line, in key lock menu item no SV)

**PV:** present value, measured by the device (lower line, in key lock menu item also upper line)

### Main switch

- On / off switch.

### Timer lamp

- Is blinking in menu items for timer modulation
- Shines continuous in the key lock menu item.

**Heating lamp**

- Shines if the heater is activated.

**Refrigorator lamp**

- Shines continuous in the key lock menu item when refrigerator runs.

**Jog-Shuttle knob**

- Change SV by turning.
- Push the knob for start and stop or probate set values.

**Safety switch**

- The safety switch turns the heating off if your set maximal temperature is overshoot.
- Set the maximal temperature rough 20 – 30 °C over the desired maximal temperature, otherwise the device maybe not reaches your set value.

**MAIN button**

- For changing between the different menu items.


**SUB button**

- For reaching the sub menu Items.

**6 Job description**

**6.1 Menu structure**

The **Wisdom WIM Multi-Room Incubator** proffer you the genial combination of MODE and SUB button and the Jog-Shuttle knob. That allows a simple and comfortable usage. The menu structure is well-arranged and pictured as follows.

**6.2 Set turn-off delay**

The turn-off delay timer serves as automatically stop function. It turns off the heating after elapsing of the timeframe you set.

- First menu item after switching the device on is the menu “set turn-off delay”. The timer

lamp is blinking.

- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- "00:00" stands for endurance run.
- Turn the Jog-Shuttle to choose your desired value and probate it by pushing the Jog-Shuttle. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Shuttle once again.

### 6.3 Set turn-on delay

The turn-on delay timer serves as automatically start function. It turns on the heating and turn-off delay timer after the timeframe you set is elapsed.

- To reach this sub menu item you have to push the SUB button in the menu "set turn-off delay"
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- Turn the Jog-Shuttle to choose your desired value and probate it by pushing the Jog-Shuttle. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Shuttle once again.
- Turn the Shuttle knob clockwise to increase the time and counter-clockwise to reduce it.
- You only reach this menu item if there is no turn-off delay time running.

### 6.4 Set temperature

- To reach this menu item you have to push the MAIN button in the menu "set turn-off delay"
- The temperature is adjustable from 0 °C to 60 °C.
- Turn the Jog-Shuttle to choose your desired value and probate it by pushing the Jog-Shuttle. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Shuttle once again.
- Turn the Shuttle knob clockwise to increase the temperature and counter-clockwise to reduce it.
- The heating lamp shines, if the heater is activated.

### 6.5 Set Offset

If the temperature measured by you differs from the temperature the device is showing (PV) you have the possibility to compensate the divergence by the Offset. So you can balance differences of about  $\pm 10^{\circ}\text{C}$ . The device will show the same temperature as the thermometer you used.

The Offset is calculated as follows:

**Temperature measured by you – showed Temperature (PV) = Offset**

- To reach this menu you have to push the SUB button in the menu item "set temperature".
- You can recognize this menu at the PV showed in the SV-line.
- If you turn the Jog-Shuttle knob.
- Set the calculated value by turning the Jog-Shuttle. The Offset will be added to the PV immediately.
- Save and confirm the Offset by pushing the MODE button, you will also reach the previous menu.

## 6.6 Key lock

- To avoid accidental changing of the set values, push the MODE button in the menu item "set temperature".
- In this menu item can you see the actually temperature in the SV-line and the remaining time of turn-off delay in the PV-line. The function of the Jog-Shuttle knob is deactivated.
- Press the MODE button once more to leave the key lock. Then you are in the menu item "set turn-off delay".

## 6.7 Storage function

The **WIM Multi-Room Incubator** includes a storage function for temperature and time values. That is practical if you want to work with the same parameters often. Then you didn't need to set them in every usage.

- For storing values you only need to hold the Jog-Shuttle knob pushed for two seconds. Successful storing is displayed by blinking of the SV for three times.
- Any stored parameters are erasable. Turn the value with the Jug-Shuttle till "0" and then hold the knob pressed for two seconds.
- If you want to overwrite a value it is not necessary to erase the old one. Set the desired value and use the Jog-Shuttle knob for storing.

## 7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.

## 8 Troubleshooting Guide

Situation	Confirmations and Solutions
Power Failure	<ul style="list-style-type: none"> <li>✓ Check the supplied voltage is present.</li> <li>✓ Check the power connection and the power cord for proper connection and integrity.</li> <li>✓ Check if the MAIN Power Switch is "ON".</li> <li>✓ Fuse is blown.</li> </ul>
Unit Not Heating	<ul style="list-style-type: none"> <li>✓ Check if the MAIN Power Switch is "ON".</li> <li>✓ Check if the set temperature (SV) confirmed with the Jog-Shuttle knob (by turn the Jog- Shuttle knob cannot change the set temperature – see point 6.4)</li> <li>✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV).</li> </ul>
Heating, but unit not reaching temperature or temperature is fluctuated	<ul style="list-style-type: none"> <li>✓ Check the Set Temperature is higher than ambient temperature (recommended min set temp: 10°C higher than ambient temp).</li> <li>✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV).</li> <li>✓ Check the Ventilation openings are closed.</li> <li>✓ Problem in the sensor or main controller.</li> </ul>
Unit Not Refrigerating (WIM-R4/ -RL4)	<ul style="list-style-type: none"> <li>✓ Check if the MAIN Power Switch is "ON".</li> <li>✓ Check if the set temperature (SV) confirmed with the Jog-Shuttle knob (by turn the Jog-Shuttle knob cannot change the set temperature – see point 6.4).</li> <li>✓ Check that the set temperature (SV) is lower than 35 °C.</li> </ul>
Refrigerating but it doesn't fall to the SV temperature (WIM-R4/ -RL4)	<ul style="list-style-type: none"> <li>✓ Problem with the the cooling circuit.</li> </ul>
ERROR 1	<ul style="list-style-type: none"> <li>✓ Temperature Sensor Error.</li> </ul>
ERROR 3	<ul style="list-style-type: none"> <li>✓ Door open alarm and appears when you open the door.</li> <li>✓ If Error 3 develops during door closed, check the Limit Switch is pressed properly.</li> </ul>

If other problems arise, or one of the above not solves problems, please contact your official agent or the manufacturer.

Reparations and amendments should only made by **witeg** authorized people. Contraventions exclude the warranty.

## 9 Technical data

		WIM-4	WIM-R4	WIM-RL4
Dimension (wxdxh)	Intern	500 x 500 x 500 mm (x 4 Rooms)		
	Extern	1.655 x 845 x 1.710 mm		
Capacity		125 l x 4 Rooms		
Heating Power		500 W x 4ea	650 W x 4ea	
Refrigerating Power		none	190 W (x 4 Rooms)	
Illumination		none		FRL 55W x 2 (x 4 Rooms)
Temp.	Range	Ambient +5°C-70°C	0 °C to 60 °C	
	Uniformity	±0,3°C at 37°C ±0,7°C at 50°C	±0,5°C at 37°C	
	Accuracy	±0,2°C at 37°C		
	Sensor	PT100		
Heat-up Time		15 min. to 37°C 30 min. to 50°C	20 Min. to 37°C 30 Min. to 50°C	
Recovery Time (Door open 30sec)		8 min. to 37°C 10 min. to 50°C	10 min. to 37°C 12 min. to 50°C	
Controller		Digital Fuzzy Controller with Jog-Shuttle Switch (Turn + Push)		
Display		Digital LCD Display with Back-Light Function		
Timer		99 hr 59 min, (Delay & Continuous Function Included)		
Material	Internal	Stainless Steel (#304)		
	External	Powder Coated Steel		
Shelves		PTFE-coated steel shelves		
Inner Door		Tempered Glass Door		
Permitted Total Load		Total Load 128 kg (16kg x 8 Shelves)		
Circulation		Forced Convection-Blower FAN		
Safety Device		Over temp. and over current protector, sensor error detector		
Others		Storage function (Temp., Timer), Locking mode (Jog-Shuttle input disabled) Alarm (Error status and timer-end)		
Net Weight		343 kg	350 kg	
Power Supply		AC 230 V, 1N~ 50/60 Hz	AC 400 V, 3~, 50/60 Hz	

## 10 Certificate

# Certificate

for

**The quality and all features were checked by the manufacturer before the shipment.**

We grant from date of purchase

**two years guaranty.**

This certificate excludes damages by natural disasters or incorrect usages by the customer.

Please look on your account and complete following table:

<b>Article</b>	<b>WIM Multi-Room Incubator</b>
<b>Typ</b>	
<b>Serialno.</b>	
<b>Date</b>	

**witeg Labortechnik GmbH**

Am Bildacker 16  
D-97877 Wertheim  
TEL: +49-9342-9301-0  
FAX: +49-9342-9301-77  
Email: info@witeg.de  
[www.witeg.de](http://www.witeg.de)