WOF

Umluftofen

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Horizontal Flow Oven

Operating Manual

gilt für / written for
WOF -50/ -105/ -155

witeg Labortechnik GmbH

(English manual at page 11)

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Contents

1 Introduction	12
2 Safety Instructions	12
2.1 Description of symbols.....	12
2.2 Instructions	12
3 Product description	14
3.1 Key features.....	14
3.2 General survey	14
3.3 Operating surface.....	15
3.4 Term description	15
4 Positioning the device and commissioning.....	16
5 Job description	16
5.1 Menu structure	16
5.2 Set turn-off delay.....	16
5.3 Set turn-on delay.....	17
5.4 Set temperature.....	17
5.5 Temperature Offset (Calibration).....	17
5.6 Key lock	17
5.7 Storage function.....	18
6 Packing contents.....	18
7 Maintenance and cleaning.....	18
8 Troubleshooting Guide	19
9 Technical data.....	20
10 Certificate.....	21

1 Introduction

Thank you for purchasing our **WOF Horizontal Flow Oven**. You did a good choice buying a **Wisdom** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a save place for future reference and hand it out if you consign the device to a third person.

Witeg reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

2 Safety Instructions

2.1 Description of symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

**ATTENTION to high temperatures**

Instructions to avoid injuries caused by high temperatures.

2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.

- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.
- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- Pay attention if you open the door, hot air is coming out.
- Do not put heat-sensitive items at the device.
- Warn all attendees if you leave the device unattended.
- Only use adequate sized vessels and avoid an over boil of liquids.
- Do not clean the scorching device.

- Reparations and amendments should only made by **witeg** authorized people. Contraventions exclude the warranty.
- Make sure, that the safety switch is adjusted to 20-30°C over max. Temperature. Otherwise the heating maybe doesn't work.

- Use this device only if you read the whole safety instructions.
- This device should only used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reach the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- All circuits and electrotechnical components are **witeg** owned technology. In fact of this it is necessary that only **witeg** authorized personal repair the device.
- Pull out the power cable if you didn't want to use the device for a longer time.

- It is important that the device stands up on a horizontal and firm subsurface to avoid vibrations.
- If you set the device in a site with vibrations, make sure that the device can not slip away.
- The device could get damaged if the vibrations are to strong.
- To avoid heat congestion or fires make sure that there is enough space for air circulations.

- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superiorly to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- Do not use this device in area with potentially explosive atmosphere.

3 Product description


3.1 Key features

The **Wisdom WOF Horizontal Flow Oven** is suitable for drying, baking, conditioning, curing, pre-heating and aging applications up to 250 °C. The device was developed for simplest usage and a very long life.


Key features are as follows:

- Simplest control by Jog-Dial knob.
- Fast commissioning.
- Long life by simply and solid construction.
- Time and temperature showed on digital LCD with Back-Light.
- Storage functions for adjusted time and temperature values.
- Turn-off and turn-on delay from one minute up to over 99 hours in one minute steps adjustable.
- Signal sound after end of delay time.
- Key lock against inadvertent changing of set temperature and time values.
- Two height-adjustable steel shelves with chrome coating included.
- Over temperature an over current protection.
- Inner room made o stainless steel.

3.2 General survey

3.3 Operating surface

3.4 Term description

SV: set value (upper line, in key lock menu items no SV)

PV: present value, measured by the device (lower line, in key lock menu item also upper line)

Main power

- On / off switch.

Timer lamp

- Is blinking in menu items for timer modulation
- Shines continuous in the key lock menu item.

Heating lamp

- Shines if the heater is activated.

Error lamp

- Is lightning if the door is opened during heating process or another error is appeared.

Jog-Dial knob

- Change SV by turning.
- Push the knob for start and stop or probate set values.

Safety switch

- The safety switch turns the heating off if your set maximal temperature is overshoot.
- Set the maximal temperature rough 20 – 30 °C over the desired maximal temperature, otherwise the device maybe not reaches your set value.

MAIN button

- For changing between the different menu items.

SUB button

- For reaching the sub menu Items.

RS232C Connector

- Monitoring and Controlling with PC by connecting RS232C cable.


4 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your customer service.
4. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.
5. Make sure that you follow all safety instructions.
6. Connect the device with the power outlet.
7. Switch the device on.
8. Set your values and start the process (look at job description).

5 Job description

5.1 Menu structure

5.2 Set turn-off delay

The turn-off delay timer serves as automatically stop function. It turns off the heating after elapsing of the timeframe you set.

- First menu item after switching the device on is the menu “set turn-off delay”. The timer lamp is blinking.
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- “00:00” stands for endurance run.
- Turn the Jog-Dial to choose your desired value and probate it by pushing the Jog-Dial. The delay time starts to run by probating the value.

- If you want to change the confirmed value, you have to push the Jog-Dial once again.

5.3 Set turn-on delay

The turn-on delay timer serves as automatically start function. It turns on the heating and turn-off delay timer after the timeframe you set is elapsed.

- To reach this sub menu item you have to push the SUB button in the menu “set turn-off delay”
- The delay time is adjustable from one minute till over 99 hours (00:01 - 99:59) in one minute steps.
- Turn the Jog-Dial to choose your desired value and probate it by pushing the Jog-Dial. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Dial once again.
- Turn the Dial knob clockwise to increase the time and counter-clockwise to reduce it.
- You only reach this menu item if there is no turn-off delay time running.

5.4 Set temperature

- To reach this menu item you have to push the MAIN button in the menu “set turn-off delay”
- The temperature is adjustable from 5 °C till 250 °C in 1 °C-steps.
- Turn the Jog-Dial to choose your desired value and probate it by pushing the Jog-Dial. The delay time starts to run by probating the value.
- If you want to change the confirmed value, you have to push the Jog-Dial once again.
- Turn the Dial knob clockwise to increase the temperature and counter-clockwise to reduce it.
- The heating lamp shines, if the heater is activated.

5.5 Temperature Offset (Calibration)

If the temperature measured by you differs from the temperature the device is showing (PV) you have the possibility to compensate the divergence by the Offset. So you can balance differences of about ± 10 °C.

The Offset is calculated as follows:

Temperature measured by you – showed temperature (PV) = Offset

- To reach this menu you have to use the SUB button in the menu item “set temperature”.
- Set the calculated value by turning the Jog-Dial. The Offset will be added to the PV immediately.
- Save the offset by holding the Jog-Dial pushed until the values in PV and SV line starts blinking.
- Push the MAIN button to leave the justification menu.

5.6 Key lock

- To avoid accidental changing of the set values, push the MAIN button in the menu item “set temperature”.
- In this menu item can you see the actually temperature in the SV-line and the remaining time of turn-off delay in the PV-line. The function of the Jog-Dial knob is deactivated.
- Press the MAIN button once more to leave the key lock. Than you are in the menu item “set turn-off delay”.

5.7 Storage function

The **WOF Horizontal Flow Oven** includes a storage function for temperature and time values. That is practical if you want to work with the same parameters often. Then you didn't need to set them in every usage.

- For storing values you only need to hold the Jog-Dial knob pushed for two seconds. Successful storing is displayed by blinking of the SV for three times.
- Any stored parameters are erasable. Turn the value with the Jog-Dial till "0" and then hold the knob pressed for two seconds.
- If you want to overwrite a value it is not necessary to erase the old one. Set the desired value and use the Jog-Dial knob for storing.

6 Packing contents

quantum	denomination
1	WOF Standard Oven
2	height-adjustable shelve
4	shelve holder
1	operating manual

7 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.

8 Troubleshooting Guide

Situation	Confirmations and Solutions
Power Failure	<ul style="list-style-type: none"> ✓ Check the supplied voltage is present. ✓ Check the power connection and the power cord for proper connection and integrity. ✓ Check if the MAIN Power Switch is "ON". ✓ Fuse is blown.
Unit Not Heating	<ul style="list-style-type: none"> ✓ Check if the MAIN Power Switch is "ON". ✓ Check if the set temperature (SV) confirmed with the Jog-Dial knob (by turn the Jog-Dial knob cannot change the set temperature – see point 5.4) ✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV).
Heating, but unit not reaching temperature or temperature is fluctuated	<ul style="list-style-type: none"> ✓ Check the Set Temperature is higher than ambient temperature (recommended min set temp: 10°C higher than ambient temp). ✓ Check the Safety switch is set higher (min. 10°C-20°C) than the set temperature (SV). ✓ Check the Ventilation openings are closed. ✓ Problem in the sensor or main controller.
ERROR 1	<ul style="list-style-type: none"> ✓ Temperature Sensor Error.
ERROR 3	<ul style="list-style-type: none"> ✓ Door open alarm and appears when you open the door. ✓ If Error 3 develops during door closed, check the Limit Switch is pressed properly.

If other problems arise, or one of the above not solves problems, please contact your official agent or the manufacturer.

Repairs and amendments should only be made by [witeg](#) authorized people. Contraventions exclude the warranty.

9 Technical data

		WOF-050	WOF-105 WOF-W105	WOF-155 WOF-W155
Dimension (wxdxh)	Internal (mm)	350x370x420	480x485x535	475x550x600
	External (mm)	518x646x722	653x711x879	718x776x944
Capacity		50 Lit.	105 Lit.	155 Lit.
Heater		650 W	1.4 KW	1.6 kW
Temp.	Range	Ambient +5 °C up to 250°C		
	Accuracy	±0.3 °C at 100°C, ±0.5 °C at 150°C		
	Uniformity	±1°C at 100°C, ±1,8°C at 150°C		
Heat-up Time		25min. to 100°C 43min. to 150°C	25min. to 100°C 40min. to 150°C	
Recovery Time (Door open 30sec)		8min. to 100°C 8min. to 150°C	6min. to 100°C 6min. to 150°C	
Controller		Digital Fuzzy Controller with Jog-Dial knob (Turn + Push)		
Display		Digital LCD Display with Back-Light Function		
Timer		99 h 59 min, (Delay & Continuous Function Included)		
Material	Internal	Stainless Steel (#304)		
	External	Powder Coated Steel		
	Insulation	Glass Wool		
Shelves		2 Stainless Steel Shelves of Adjustable Position		
Door Gasket		High Temperature Grade Foamed Silicone Rubber		
Viewing Window		Tempered Safety Glass, 3 Fold, 5t		
Ventilation		Stainless Steel Cap, Dia. 40mm		
		1ea	2ea	
Sensor		PT100		
Circulation		Forced Convection-Blower FAN		
Safety Device		Over Temp. and Over Current Protector, Sensor Error Detector		
Others		Storage Function (Temp., Timer), Locking MAIN (Jog-Dial Input disabled) Alarm (Error Status and Timer-end)		
Net Weight (kg)		44	69	78
Packing Size (wxdxh)		622x754x807 mm	743x819x1.068 mm	808x884x1.133 mm
Shipping Weight(kg)		59	79	92
Power Consumption		674 W	1390 W	1566 W
Power Consumption	at 100°C	117 Wh	187 Wh	204 Wh
	at 150°C	230 Wh	337 Wh	398 Wh
Power Supply		AC 230 V, 1N~, 50/60 Hz		

10 Certificate

Certificate

for

The quality and all features were checked by the manufacturer before the shipment.

We grant from date of purchase

two years guaranty.

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	WOF Horizontal Flow Oven
Typ	
Serialno.	
Date	

witeg Labortechnik GmbH
Am Bildacker 16
D-97877 Wertheim
TEL: +49-9342-9301-0
FAX: +49-9342-9301-77
Email: info@witeg.de
www.witeg.de