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# 1 Introduction

Thank you for purchasing our **HT-T Overhead Stirrer with Remote controller**. You did a good choice buying a **Wisdom** product.

Before use, please read this manual and use it properly.

To cover the user's health and avoid the product's destruction you have to follow the safety instructions of the second chapter.

Please keep this operating manual in a safe place for future reference and hand it out if you consign the device to a third person.

**Witeg** reserves to modify the capsule or make amendments for the improvement of quality and efficiency without prior notice.

## 2 Safety Instructions

### 2.1 Description of Symbols

**WARNING**

Information for avoiding injuries or fatal accidents.

**ATTENTION**

Safety instruction to avoid a damage of the product.

**ATTENTION**

Information about technical parameters.

**Environment**

Information about the operating environment.

**WARNING of explosions**

Safety instruction to avoid explosions.

### 2.2 Instructions

For a safety and faultless using of the device it is indispensable that you read the following instructions purposely and follow them by the use.

- Connect the device only to the suitable power supply. Please have a look to the identification plate.
- Use the device only in electric circuits with ground fault circuit interrupter.
- Check the power supply cord before you stick it into the electrical outlet.
- Do not use power bars to extend the power supply cord.

- Only use extension cables with sufficient cable cross section.
- Do not touch the device with wet hands during the usage because this might cause an electrical shock.
- If unexpected sound, smell or smog is generated by the device, pull out the main plug and contact the manufacturer or your supplier immediately.
- Disconnect the power supply plug before cleaning or maintenance.
- If you use radiation or contaminated samples the responsibility is only up to you.

- Use this device only if you read the whole safety instructions.
- This device should only be used by persons who are familiarized with safety instructions made for laboratories.
- Pay attention that no liquid reaches the controller or inward of the capsule.
- You have to keep the device clean if you want to get a longer life of it.
- Wear corresponding protective work clothing during the usage.
- Pull out the power cable if you didn't want to use the device for a longer time.

- Do not use this device in areas with potentially explosive atmosphere.

- Make sure that the device can not slip away.
- The device could get damaged if the vibrations are too strong.
- To avoid heat congestions or fires make sure that there is enough space for air circulations.
- Avoid direct solar radiation.
- Do not use the device near highly combustible materials.
- Do not set the device outdoors.
- It is superior to set the device in a good illuminated place to avoid accidents engendered by handling errors.
- To set the device in wet or dusty places could cause overheating, short circuits or fires.

- All circuits and electro technical components are **witeg** owned technology. In fact of this it is necessary that only **witeg** authorized personal repair the device. Contraventions exclude the warranty.

## 3 Product description


### 3.1 Key features

The **HT-T Overhead Stirrer with Remote Controller** is designed for chemical, biological, medical and pharmaceutical applications. With the possibility of remote controlling the device is designed for usage under a fume hood or behind a splash-guard. For controlling the device no direct contact to the device is needed.

Key features are as follows:

- Long life by simply and solid construction.
- Controlling with remote controller
- Jog-Shuttle knob for easy and fast handling.
- Up to 1.000 rpm.
- Max. viscosity: 150.000 mPas.
- Constant speed in spite of viscosity fluctuation.
- Motor provided with overload and overheat protection.
- Stepless variable adjustability of every stirrer concerning to the continuous bored chuck.
- Power bar graph and speed showed on illuminated LCD.
- Smooth increase of speed for higher safety.
- Storage functions for speed value.

### 3.2 General survey

### 4 Packing contents

Units Only			Packages	
quantity	DH.WOS01065	DH.WOS01075	DH.WOS01066	DH.WOS01076
1	HT-50T	HT-120T	HT-50T	HT-120T
1	Allen key			
1	body holder			
1	chuck key			
1	operating manual			
1			clamp <b>CL200</b>	
1			impeller <b>PL020</b>	
1			stand rod <b>ST120</b>	stand rod <b>ST200</b>

## 5 Positioning the device and commissioning

1. Follow especially the safety instructions marked with following symbol:
2. Remove the package and perhaps soiling carefully.
3. If there are any unexpected transport damages, please inform your forwarder and your customer service.
4. Attach the clamp to the tripod rod.
5. Stick the device at the clamp and tighten the screw.
6. Check that each connection fits close.
7. Fix the stirrer in the chuck with included chuck key.
8. connect the remote controller to the device
9. Before you connect the device with the electrical outlet, please let some time laps until the device is acclimatised.
10. Make sure that you follow all safety instructions.
11. Connect the remote controller with the power outlet.
12. Switch the device on and start the process.

## 6 Job description

- First switch the device on.
- If there is a value stored the stirring process starts by pushing the Jog-Shuttle.
- If there is no value stored the stirring process starts by turning the Jog-Shuttle.
- Change the set speed during the stirring application by turning the Jog-Shuttle.
- Start and stop by pushing the Jog-Shuttle.

## 7 Storage function

- If you work with the same speed often you can store this value and thereby you need not to set often used speed values in every usage.
- To store values you only have to hold the Jog-Shuttle knob pushed for two seconds until the value in the first line in display is blinking.
- Delete stored values by over-storing with "0".

## 8 Maintenance and cleaning

- Please read the whole safety instructions before you start to maintain or clean the device.
- Pull the power plug out of the electrical outlet.
- Only use care cleaning agents and soft cleaning rags.
- Organic lotions, strong chemicals and rough cleaning rags can damage the device.
- If you not want to use the device for a long time, pull of the power plug and keep it dry in the package.
- Only use original spare parts.


## 9 Technical data

		HT-50T	HT-120T
<b>Stirring Capacity (H<sub>2</sub>O)</b>		Max. 40 Lit.	Max.60 Lit.
<b>Max. Viscosity (mPas)</b>		100.000	150.000
<b>Motor Rating (Input / Output)</b>		83 / 75	130 / 120 W
<b>Maximum Torque (Ncm)</b>		300	500
<b>Max. Shaft Diameter (mm)</b>		10	
<b>Speed Range (rpm)</b>		50 – 1.000	
<b>Controller</b>		Digital Feedback Controller with Jog Shuttle knob (Turn + Push)	
<b>Display</b>		Digital LCD Display (RPM and Power bar-graph)	
<b>Operating Environment</b>	<b>Ambient Temp. Range (°C)</b>	0 ~ 40	
	<b>Max. Permissible Humidity (%)</b>	80	
<b>Safety Device</b>		Over Load & Heat Protector for motor	
<b>Dimension (wxdxh mm)</b>		155x180x110	155x230x110
<b>Device-Weight (kg)</b>		3,5	4
<b>Others</b>		<ul style="list-style-type: none"> <li>- Smooth Start and Change of the Stirring Speed</li> <li>- Storage Function (RPM)</li> <li>Alarm (Error status)</li> </ul>	
<b>Power Supply</b>		AC 230 V, 1N~, 50/60 Hz	

## 10 Certificate

# Certificate

for

**The quality and all features were checked by the manufacturer before the shipment.**

We grant from date of purchase

**two years guaranty.**

This certificate excludes damages by natural disasters or incorrect usages by the costumer.

Please look on your account and complete following table:

Article	HT-T Overhead Stirrer
Typ	
Serialno.	
Date	

**witeg Labortechnik GmbH**

Am Bildacker 16  
D-97877 Wertheim  
TEL: +49-9342-9301-0  
FAX: +49-9342-9301-77  
Email: info@witeg.de  
[www.witeg.de](http://www.witeg.de)