



K905XX **AUTOMATIC POTENTIOMETRIC TITRATOR**

OPERATION AND INSTRUCTION MANUAL

Koehler Instrument Company, Inc.

1595 Sycamore Avenue • Bohemia, New York 11716-1796 • USA

Toll Free: 1-800-878-9070 (US only) • Tel: +1 631 589 3800 • Fax: +1 631 589 3815

<http://www.koehlerinstrument.com> • e-mail: info@koehlerinstrument.com

Petroleum Testing & Analysis Instrumentation • Custom Design & Manufacturing

CERTIFICATE OF CONFORMANCE


Automatic Potentiometric Titrator K905XX

This certificate verifies that part number K905XX, Automatic Potentiometric Titrator, was manufactured in conformance with the applicable standards set forth in this certification.

Specifications:

ASTM D664
ASTM D2896
ASTM D3227
ASTM D4739

This unit is tested before it leaves the factory, to ensure total functionality and compliance to the above specifications and ASTM standards. Test and inspection records are on file for verification.



Jesse Kelly
Application Engineer
Koehler Instrument Company

CONTENTS

SECTION 1

Introduction	5
Specifications	5
Features	6

SECTION 2

Familiarization with instrument	9
Front Panel – Keyboard Layout and Functions	10

SECTION 3

Instrument Operation - Switching 'ON' the instrument	12
1) Develop Method	
♦ Generating the Quick Method	13
♦ Viewing & Running & Copying the Last Quick Method	15
♦ Generating the Full Method	19
♦ Viewing, Copying, Skipping & Running Full Method	41
2) Function	
♦ Standardization & Calibration	45
♦ mV indicator & Priming	47
♦ Priming	47
♦ Burette Replacement	48
♦ Manual Titration	49
3) Configuration	
♦ Setting RTC (Clock)	50
♦ Balance Connectivity	50
♦ Temperature Compensation	50
♦ Printer Selection	51
♦ Password Entry	51
♦ Data Transfer	52
♦ Resolution Selection	53
♦ Titrant Temp Factor	53
4) Print	
♦ Printing Last Titration Result	54
♦ Printing Data Table	55
♦ Printing The Titration Analysis Condensed Report	55
♦ Printing the Titration Statistics Report	55
♦ Printing pH Calibration Report	55
♦ Printing Method Parameters	56
♦ Printing mV v/s μ l Grapha Report	56

5) Display		
♦ Display Result	57
♦ Display Statistics	57
♦ Display Calibration Data	57
6) Reprocess		58
7) Run		
♦ Run Last Method	59
♦ Run New Method	60
♦ Run Last Quick Method	61
SECTION 4		
Result Calculation	64
Entering Blank Volume/Primary Dose	65
Entering Constant related to Sample	66
Automatic determining and storing of molarity	68
Maintenance and cleaning of burette	68
Maintenance and cleaning of electrode	69
Appendix A : Alphanumeric entry	70
Appendix B : Error Codes	71
Appendix C : Printer DIP switch setting	73
Appendix D : Parameter selection table	74
Appendix E : Titration Type	76
Appendix F : Titration Modes	77
Appendix G : Method Parameters	78
Appendix H : Sample Analysis Parameters	79
Appendix I : Calculation Parameters	81
Appendix J : Standard Methods (41-50)	82
Appendix K : PC Setup instruction for Data Downloading	87

.....

.....

THE CONTENTS OF THIS MANUAL ARE SUBJECT TO CHANGE WITHOUT NOTICE.

SECTION 1

- **Introduction**
- **Specifications**
- **Features**

K905XX- AUTOMATIC POTENTIOMETRIC TITRATOR

INTRODUCTION

An Automatic Titrator is a combination of electrochemical potential measurement through electrode with microprocessor driven burette, thus increasing the accuracy of measurement. The instrument is capable of performing almost all basic titrimetric functions and cater to all types such as 1) Aqueous 2) Non-aqueous 3) Redox 4) Precipitation / Argentometric 5) Complexometric / EDTA 6) Voltametric 7) Back titration with simplicity. User can define particular method for a specific titration.

SPECIFICATIONS :

Control	: Micro-controller based.
Burette	: 1ml, 5 ml or 10 ml capacity, interchangeable. 25 ml Optional.
Burette resolution	: 1/10000 for 10 ml , 1/5000 for 5 ml and 1/1000 for 1 ml.
Filling time	: < 20 seconds.
Principle	: Volume determination by equivalence point, end point or pH STAT.
Modes of operation	: a)Incremental b)Equilibrium c)Cut-off by pH d)pH STAT titration.
Methods	: a)Acid-base b)Non-aqueous c)Redox d)Precipitation e)Complexometric f)Back g)Photometric h) Thermometric.
mV range	: ± 3200 mV
Accuracy	: ± 0.1 mV (± 0.0016 pH).
Amplifier Input Impedance	: $> 10^{12}$ ohms.
End Point detection	: a)Potentiometric b)Voltametric c)Thermometric and d)Photometric.
Cut-off criteria	: a) Volume b)End point c)mV/pH.
Results	: Molarity, % Assay (wt), % volume (ml) , ppm , mg/l, mg/g, g/l, meq/l, mol/kg, TAN and TBN for oil samples.
Sensors	: 1)Electrodes for Potentiometric titration - pH, Ion, Redox, Precipitation. a)Any combination electrode with BNC connector, b)Differential System comprising Sensing (Indicator) with BNC and Reference (4mm banana connector). 2) Electrodes for KF/Voltametric titration with TNC connector. 3) Temperature (PRT/PT100) for ATC with 5-pin shell connector
Calibration	: 3-point Calibration with user entered buffer values and standardization (7pH)
Titration Head	: Manual stand with swiveling arm for electrode, dispensing tip & stirrer is provided as standards
Stirring System	: Micro-controller based variable speed, high torque stirrer with indication.
Data Storage	: Non-volatile memory.
Method Storage	: 50 methods with parameters.
Titrant Molarity storage	: within Method data.
Input/Output Peripheral Interface	: a) Parallel Port : 1 No. for printer. b) Serial Port : 2 Nos. for Balance & PC.
Report Format	: Selection of various formats - a) Method parameters b) Titration analysis report c) Titration analysis condensed report d) Titration data table e) Titration graphic report - 1) $\mu\text{l v/s mV}$ 2) $\mu\text{l v/s First derivative}$ 3) $\mu\text{l v/s Second derivative}$ 4) $\mu\text{l v/s time}$. f) Statistics report g) Calibration report.
Keyboard	: Alphanumeric splash water proof polyester soft keys.
Display	: 40 x 2 line back lighted LCD display.
Convertibility	: Can be used KF titrator for Moisture estimation with KF burette system (Optional)

Environmental Operating Conditions :

- a) Operation : Indoor.
- b) Temperature : Ambient to 45 °C
- c) Operating Temperature and humidity : +5 to 45°C Max, 5 to 90% non condensing.
- d) Storage Temperature and humidity : +5 to 45°C Max, 5 to 90% non condensing.
- e) Altitude in meters : up to 2000m

Power : 110V/220V AC ±10%, 60/50Hz. 40W.

Weight and dimensions :

- Control Unit : 12 inch x 12 inch x 6 inch (W X D X H)
- Titration unit : 12 x 16 x 16 (W X D X H)
with Vortex Stirrer
- Magnetic Stirrer : 6 x 10 x 10 (W X D X H)
(For KF Titration - Optional)
- System Weight : 9 Kg.

FEATURES:

- Advanced Micro-controller based user-friendly state-of-the-art product design with alphanumeric splash waterproof polyester soft keys for keyboard and User interactive software in dialogue mode for ease of operation.
- Quick interchangeable burette assemblies with intelligent recognition for its volume size. Burette factor for dispensing correction is available for true end point volume.
- Alphanumeric entry of Sample Name, Titrant Name, Identification Number, Date & Time with type of electrode used for authentication. Daily Auto Incremented Run number and Factory entered CUSTOMER NAME & Instrument Sr.No. on report printouts make the system foolproof and GLP compliant.
- Protection against invalid parameter entries during method generation.
- Storage of 40 User methods with parameters. (Plus10 not editable – Default methods)
- Two Tier - ADMIN & USER Password Protection for Method data during ADD/EDIT, DELETE, COPY, and for Calibration.
- Three modes of titration - incremental, equilibrium and cut-off by pH mode to perform almost all types of titration. By selecting titration method, instrument prints the type of appropriate electrode.
- Automatic selection of Sensor Input as per titration method selected.
- During titration, the measured variable i.e. electrode potential (mV) or pH value is shown on display together with dispensed volume and number of End Point (EP) detected.
- User selectable End Point (EP) evaluation up to 9 EP during the run, and calculation by first, last, largest, all or selected EP with display of results and printout.
- Composite Differential Electrode Amplifier unit for potentiometric and voltametric / KF titration, having BNC connectivity to various electrodes and temperature sensor for indication.
- Capability to use as simple pH meter with display of pH, mV values & Temperature display.
- Facility to use as a dispenser for fixed volume dosing or dilution allows to perform manual titration with user defined dose and mV indication.
- User Programmable selectivity for report format, complying with GLP requirements: a) Report giving titration parameter and result. b) Data table giving mV, pH, δ mV, mV/ml, 2nd Deriv. and volume. c) Graphics report giving mV v/s μ l titration curve. d) Graphics report 1st deriv. graph v/s μ l titration curve e) Graphics report of 2nd deriv curve. f) Report of method parameters for 50 methods. g) Condensed report of titration parameters and results. h) Auto evaluation report for multi EP samples - EP1, EP2-EP1 & EP3-EP2 etc. available. The reports can be obtained even after resetting/power off / power fail conditions or Run Aborted.

FEATURES:**contd...**

- Result calculation facility to obtain printout in different units such as molarity, factor, % assay (wt), % volume (ml), ppm, mg/l, mg/g, ml/g, g/l, mEq/l, mol/kg, TAN & TBN.
- Reprocessing of threshold and recalculation of EP without performing the new run.
- Display and printout of daily Auto Incremented Run (AIR) number which makes system foolproof and GLP compliant.
- Statistic function with run selectivity for finding Mean, SD, RSD & CV of last 10 repeat run results could be viewed or printed.
- Real Time Clock (RTC) for time display and report printout with run time indication.
- Special Method Compliant to ASTM D664, D2896 & D4739 for TAN and TBN analysis for oil samples.
- Automatic evaluation of molarity determination of titrant.
- Quick Method for Trial analysis RUN having only titration parameter.
- Error indication on Display, helps user to trace the problem.
- Titration run can be started with last run parameters.
- Display of last run or default parameters while entering or running method.
- Display of stir time after pre-dispensing (with selectable dose and time).
- Escape and Back keys for better user interaction.
- Vortex stirrer for homogeneous stirring with glass propeller (provides chemical inertness).
- Digital Speed control of vortex and magnetic stirrer with indication.
- Availability of optional in situ accessory for Titration Vessel heating to perform titration from ambient+10°C to 80 °C. Cooling can be achieved with external chiller.
- PC Compatibility – RS232 interface for Data Downloading to Result Data Capturing Software (Ethernet based) or Windows Hyper Terminal Software (Optional).
- Balance interface for automatic sample weight transfer.
- K905XX can be converted to perform Karl Fischer titration by just changing burette assembly.

SECTION 2


- **Familiarization with Instrument**
- **Front Panel – Keyboard Layout and Functions**

K905XX - Automatic Potentiometric Titrator

Familiarization with the Instrument : -

- 1) Titration Controller: System Control unit for controlling the Titration Unit.
- 2) Titration Unit : Consists of Burette Assembly having Gas tight syringe with teflon tipped plunger with luer fitting. The miniature 3-way valve is used to control the titrant flow. Two Teflon tubings are fitted with finger tight fittings. The left side tubing is for titrant aspiration and the right side tubing is for dispensing the titrant in reaction vessel.
- 3) Stirrer Assembly : Consists of Vortex stirrer with glass propeller with digital controller for variable speed regulation.
- 4) Printer Output : 25 Pin D type connector – Centronic Parallel port interface
- 5) Display : 40 X 2, Back lighted, LCD display.
- 6) Key Board : Multi-function tactile keyboard as shown below.
- 7) Electrode Connection: Differential input with BNC and 4mm Banana Connector.
- 8) Temperature sensor : 5 pin Shell connector
- 9) Balance connectivity : RS232C – DB9.
- 10) PC Connectivity : RS232C – DB9 for data downloading to Result Data Capturing Software (Ethernet based) (Optional) or Windows HyperTerminal or any other terminal software.

FRONT PANEL KEYBOARD LAYOUT :










Key functions:

L.H.S. Keys: Alphanumeric Dual Function (A~Z, 1~9), Zero (0) with Dash (-), Dot (.), Slash (/), Decimal point (.) and Clear Entry(CE) Key.

R.H.S. Keys: Functional Keys and RESET key.

Description of Functional Keys:

- 1]  : To reset the instrument or re-initialize before RUN.
- 2]  : To go back to previous entry.
- 3]  : To Return to Main menu from any operation and to abort the titration during run.
- 4]  : a) To acknowledge functions.
b) To enter titration parameters.
c) To stop continuous monitoring of mV read by electrode in Scan mode.
- 5]  : To cancel the wrong entries.
- 6]  : To enter decimal point during numeric entries. To Add / Remove the Beep during run.
- 7]  : To Start / Stop Stirrer Manually.

SECTION 3

- **Instrument Operation**
- **Titration Method Development**
- **Function**
- **Methods**

OPERATION :

- **Switching 'ON' the instrument and connecting the accessories:**

- Connect interface cable between Controller Unit & Titration Unit.
- Connect the DC Adapter power cord to 230V AC/50 Hz, 5A socket of stabilized power supply with proper earth connection & insert the DC jack in the Controller Unit socket.
- Connect pH Electrode on the back side of the instrument.
- Connect printer through appropriate cable to the backside of the instrument.
- Connect the Vortex stirrer to the backside of Stirrer control unit.
- Turn the power switch 'ON', with beep sound the LCD display shows the Model & the firmware version with :

AUTO TITRATOR Vx.x

Display

AT 10 ml Burette Assembly detected

Display

Wait...

Display

MAIN MENU

A

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

- 1) Develop Method :- To generate(add) Edit , View Titration Method parameter, Delete, Copy Method. (Refer pages 13)
- 2) Function :- To perform pH Calibration, mV scan, Burette Priming, Burette change or Manual titration. (Refer pages 45)
- 3) Config :- System Settings – Clock/date, Balance, Temperature compensation, Printer, Password setting, Titrant temperature factor., Data Transfer, Resolution (Refer pages 50)
- 4) Print :- Titration report, Statistics, Method parameter, pH Calibration report. (Refer pages 54)
- 5) Display :- Titration result, Statistics, pH Calibration data. (Refer pages 57)
- 6) Reprocess :- Data Reprocessing (Refer pages 58)
- 7) Run :- To Perform titration with desired method (Refer pages 59).

1. Titration Method Development

- Developing the Quick Method : Follow the below mentioned instructions

MEAS. TYPE	: Potentiometric	TITR. TYPE	: Acid-Base Aqueous
TITR. MODE	: Incremental	P.D.VOLUME	: 0 µl
STIR TIME	: 5 sec	INITIAL DOSE	: 100 µl
INTERVAL	: 5 sec	VOLUME LIMIT	: 10000 µl
TITR. END CRIT.	: End Point	NO. OF EP	: 1
CAL. BY	: 1 st EP	DERIV.THRESH.:	50

1)Develop Method 2) Function 3)Config. : 7
4)Print 5) Display 6) Reprocess 7)Run

1
ABC

1)Quick method 2)Full method :

1
ABC

1)Last Method 2)Edit Method : 1
3)Copy Method

2
DEF

Measurement Type
1)Potentiometric 2)Voltametric :

1
ABC

ENTER

Titration Type 1)Aq. 2)Non-Aq. 3)Redox
4)Pr. 5)EDTA 6)Back 7)Photo 8)Thermo : 1

1
ABC

ENTER

Select titration mode
1)Incr. 2)Equil. 3)Cutoff(pH) 4)pHstat: 1

1
ABC

ENTER

Predispense Volume (0-95000 µl) : 0

ENTER

Predispense Dose (5-10000 µl) : 1000

ENTER

Predispense Interval (0-999 sec) : 0

ENTER

Stir Time (0-999 sec) : 5

5
MNO

ENTER

Initial Dose (5-2000 µl) : 50

1
ABC

0
-./

0
-./

Initial Dose (5-2000 μ l) : 100 ENTER

Interval (1-999 sec) : 5 5
MNO ENTER

Volume Limit (Max. : 99999 μ l) : 1 0 0
ABC -./ -./

0 0
-./ -./ Volume Limit (Max. : 99999 μ l) : 10000 ENTER

Titration End Criteria
1] EndPoint (EP) 2] VolumeLimit (VL) : 1 1
ABC ENTER

Number of EP (1-9) : 1
ABC ENTER

Calculation by
1)1st 2)Largest 3>Last 4)All 5)Select : 1 1
ABC ENTER

Derivative Threshold (1-32000) : 5 0
MNO -./


Derivative Threshold (1-32000) : 50 ENTER

Quick Method Stored ENTER


1)Run 2)View :

- **Viewing & Running the Last Quick Method :** Follow the below mentioned steps . Details of last quick method to be viewed & run :

MEAS. TYPE	: Potentiometric	TYPE	: Acid-Base Aqueous
MODE	: Incremental	PREDISPENSE VOLUME	: 0 µl
PREDISPENSE DOSE	: 1000	PREDISPENSE INTERVAL	: 0
STIR TIME	: 5 sec	INITIAL DOSE	: 100 µl
INTERVAL	: 5 sec	VOLUME LIMIT	: 10000 µl
TITR. END CRIT.	: End Point	NO. OF EP.	: 1
CAL. BY	: 1 st EP	DERIV.THRESH.	: 50



- To repeat the same quick method with same parameters, select '1' option.
- To print the data table report for quick method, select '3' option.
- To get result for the quick method enter sample analysis parameters as following by selecting '2' option.



0	0	Conc. (0.00001-99.9999) : 0.100	ENTER
-. /	-. /		

1)BlankVolume 2)PrimaryDose 3)None :	3 GHI	ENTER
--------------------------------------	----------	-------

Constant Code (1-8) :	1 ABC	ENTER
-----------------------	----------	-------

Const 1 (0.0001-999999.99) :	6 PQR	3 GHI
------------------------------	----------	----------

Const 1 (0.0001-999999.99) : 63	ENTER
---------------------------------	-------

Sample Vol. (0.0001 – 999.99999 ml): 1.0000

Vol. Reqd : 0.9425 ml R : 0.1924 MOLARITY(S) 1)Repeat Run 2)Reprocess 3)Print :
--

- a) To repeat the same quick method with same parameters, select '1' option.
- b) To print the data table report for quick method, select '3' option.
- c) To reprocess/ recalculate result in different result unit or to evaluate endpoint by changing the threshold , select '2' option.

Derivative Threshold (1-32000) :	2 DEF	0 -. /
----------------------------------	----------	-----------

Derivative Threshold (1-32000) : 20	ENTER
-------------------------------------	-------

Calculation by 1)1 st 2)Largest 3)Last 4)All 5)Select : 1	ENTER
---	-------

Conc. (0.00001-99.9999) :	0 -. /	.	1 ABC
---------------------------	-----------	---	----------

0	0	Conc. (0.00001-99.9999) : 0.100	ENTER
-. /	-. /		

1)BlankVolume 2)PrimaryDose 3)None :	3 GHI	ENTER
--------------------------------------	----------	-------

Constant Code (1- 8) :

1
ABC

ENTER

Const 1 (0.00001-999999.99) :

6
PQR

3
GHI

Const 1 (0.00001-999999.99) : 63

ENTER

Sample Vol. (0.0001 – 999.9999 ml) : 1.0000

Vol. Reqd : 0.9425 ml R : 0.1924 MOLARITY(S)
1)Repeat Run 2)Reprocess 3)Print :

• Viewing Last Quick Method :

1)Quick method 2)Full method :

1
ABC

1)Last Method 2)Edit Method : 1
3)Copy Method

1
ABC

1)Run Method 2)View :

2
DEF

ENTER

Measurement Type :
Potentiometric

ENTER

Type : Acid-Base Aqueous
Mode : Incremental

ENTER

Predispense Volume (μ l) : 100
Predispense Dose (μ l) : 50

ENTER

Predispense Interval (sec) : 0
Stir Time (sec) : 0

ENTER

Initial Dose (μ l) : 100
Interval (sec) : 3

ENTER

VolLimit (μ l) : 10000
Titr. End Crit. : EndPoint

ENTER

No. of EP : 1 Cal. By : 1st

ENTER

Deriv. Threshold : 15

ENTER

1)Run 2)View :

• **Copying Last Quick Method :** Follow the below mentioned instructions

1)Quick method 2)Full method :

1
ABC

1)Last Method 2)Edit Method : 1
3)Copy Method

3
GHI

ENTER

Method No. (1-40) :

1
ABC

ENTER

Method No. 1 copied from Quick Method

After 4-5 seconds

1)Last Method 2)Edit Method : 1
3)Copy Method

1)Run 2)View :

- **Developing the Full Method :** Follow the below mentioned steps

Examples: TITRATION BY INCREMENTAL METHOD

METHOD NO.	: 1	TITRANT NAME	: NAOH
SAMPLE NAME	: KHP	TITR. TYPE	: Acid-Base Aqueous
MEAS. TYPE	: Potentiometric	P.D.CRITERIA	: 0 mV
TITR. MODE	: Incremental	P.D.DOSE	: 100 µl
P.D.VOLUME	: 500 µl	STIR TIME	: 5 sec
P.D.INTERVAL	: 5 sec	INTERVAL	: 5 sec
INITIAL DOSE	: 100 µl	TITR. END CRIT.	: End Point
VOLUME LIMIT	: 10000 µl	CALCULATION BY	: 1 st
NO. OF EP	: 1	BV/PD/None	: None
DERIV.THRESH.	: 50	CONSTANT	: 1
CONC. (mol/l)	: 0.100		
CONSTANT CODE	: 1		
PRINT REPORT	: None		

1)Develop Method 2) Function 3)Config. : 7 4)Print 5) Display 6) Reprocess 7)Run	1
	ABC

1)Quick method 2)Full method :	2
	DEF

1)Add/Edit method 2) Copy method 3) Skip method 4) View :	1
	ABC

Enter User Password: Enter ADMIN or USER Password (ref. Page. 51)

Method No. (1-40) :	1	ENTER
	ABC	

If method is already present , Display reads :

Method exists !! Edit 1)Y 2)N :	1
	ABC

SAMPLE NAME :	ENTER
---------------	-------

TITRANT NAME :	ENTER
----------------	-------

Enter Sample and Titrant Name [refer Appendix-A for alphanumeric entries]

Measurement Type 1)Potentiometric 2)Voltametric :	1	ENTER
	ABC	

Titration Type 1)Aq. 2)Non-Aq. 3)Redox 4)Pr. 5)EDTA 6)Back 7)Photo 8)Thermo : 1	1	ENTER
	ABC	

Select titration mode 1)Incr. 2)Equil. 3)Cutoff(pH) 4)pHstat: 2	1	ENTER
	ABC	

<

Predispense Criteria (0-1000 mV) : 0	0 -./	ENTER
--------------------------------------	----------	-------

Predispense Volume (0-95000 µl) : 0	5 MNO	0 -./	0 -./
-------------------------------------	----------	----------	----------

Predispense Volume (0-95000 µl) : 500		ENTER
---------------------------------------	--	-------

Predispense Dose (5-1000) :	1 ABC	0 -./	0 -./
-----------------------------	----------	----------	----------

Predispense Dose (5-1000) : 100		ENTER
---------------------------------	--	-------

Predispense Interval (0-999 sec) : 0	5 MNO	ENTER
--------------------------------------	----------	-------

Stir Time (0-999 sec) : 5	5 MNO	ENTER
---------------------------	----------	-------

Initial Dose (5-2000 µl) : 100	1 ABC	0 -./	0 -./
--------------------------------	----------	----------	----------

Initial Dose (5-2000 µl) : 100		ENTER
--------------------------------	--	-------

Interval (1-999 sec) : 5	5 MNO	ENTER
--------------------------	----------	-------

Volume Limit (Max. : 99999 µl) : 10000	1 ABC	0 -./	0 -./
--	----------	----------	----------

0 -./	0 -./	Volume Limit (Max. : 99999 µl) : 10000	ENTER
----------	----------	--	-------

Titration End Criteria 1] EndPoint (EP) 2] VolumeLimit (VL) : 1	1 ABC	ENTER
--	----------	-------

Number of EP (1-9) :	1 ABC	ENTER
----------------------	----------	-------

Calculation by 1)1 st 2)Largest 3>Last 4)All 5)Select : 1	1 ABC	ENTER
---	----------	-------

Derivative Threshold (1-32000) :	5 MNO	0 -./
----------------------------------	----------	----------

Derivative Threshold (1-32000) : 50		ENTER
-------------------------------------	--	-------

Conc. (0.00001-99.99999 M1) :

Conc. (0.00001-99.99999) : 0.100

1)BlankVolume 2)PrimaryDose 3)None :

Constant Code (1-20) :

Const 1 (0.00001-999999.9) :

Const 1 (0.00001-999999.9) : 1

Print Report 1)Full 2)Select 3)Condensed 4)None : 4

Method No. 1 Stored

After 4-5 seconds

1)Run 2)View :

• **TITRATION BY EQUILIBRIUM METHOD:**

METHOD NO.	: 2	TITRANT NAME	: NAOH
SAMPLE NAME	: KHP	TITR. TYPE	: Acid-Base Aqueous
MEAS. TYPE	: Potentiometric	P.D.CRITERIA	: 0 mV
TITR. MODE	: Equilibrium	P.D.DOSE	: 100 µl
P.D.VOLUME	: 500 µl	STIR TIME	: 5 sec
P.D.INTERVAL	: 5 sec	EQ. FACTOR	: 5
INITIAL DOSE	: 100 µl	TITR. END CRIT.	: Volume Limit
VOLUME LIMIT	: 10000 µl	CALCULATION BY	: 1 st
MIN.CTRL.DOSE	: 5 µl	BV/PD/None	: None
DERIV.THRESH.	: 500	CONSTANT	: 1
CONC. (mol/l)	: 0.100	T. DATA TABLE	: 2 (No)
CONSTANT CODE	: 1	DERIV. GRAPH	: 2 (No)
TITR. REPORT	: 1 (Yes)	2nd DERIV. GRAPH	: 2 (No)
mV/µl GRAPH	: 1 (Yes)		
Time/µl GRAPH	: 2 (No)		

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run	1 ABC
---	----------

1)Quick method 2)Full method :	2 DEF
--------------------------------	----------

1)Add/Edit method 2) Copy method 3) Skip method 4) View :	1 ABC
--	----------

Enter User Password:	Enter ADMIN or USER Password (ref. Page. 51)
----------------------	---

Method No. (1-40) :	2 DEF	ENTER
---------------------	----------	-------

SAMPLE NAME :	ENTER
---------------	-------

TITRANT NAME :	ENTER
----------------	-------

Enter Sample and Titrant Name [refer Appendix-A for alphanumeric entries]

Measurement Type 1)Potentiometric 2)Voltametric :	1 ABC	ENTER
--	----------	-------

Titration Type 1)Aq. 2)Non-Aq. 3)Redox 4)Pr. 5)EDTA 6)Back 7)Photo 8)Thermo : 1	1 ABC	ENTER
--	----------	-------

Select titration mode 1)Incr. 2)Equil. 3)Cutoff(pH) 4)pHstat: 2	2 DEF	ENTER
--	----------	-------

Predispense Criteria (0-1000 mV) : 0	0 -./	ENTER
--------------------------------------	----------	-------

Predispense Volume (0-95000 μ l) : 0	5 MNO	0 -./	0 -./
Predispense Volume (0-95000 μ l) : 500	ENTER		
Predispense Dose (5-1000) :	1 ABC	0 -./	0 -./
Predispense Dose (5-1000) : 100	ENTER		
Predispense Interval (0-999 sec) : 0	5 MNO	ENTER	
Stir Time (0-999 sec) : 5	5 MNO	ENTER	
Initial Dose (5-2000 μ l) : 100	1 ABC	0 -./	0 -./
Initial Dose (5-2000 μ l) : 100	ENTER		
Interval (1-999 sec) : 1	1 ABC	ENTER	
Interval (1-999 sec) : 1			
Volume Limit (Max. : 99999 μ l) :	1 ABC	0 -./	0 -./
0 -./	0 -./	Volume Limit (Max. : 99999 μ l) : 10000	ENTER
Titration End Criteria 1] EndPoint (EP) 2] VolumeLimit (VL) : 2	2 DEF	ENTER	
Calculation by 1)1 st 2)Largest 3>Last 4)All 5)Select : 1	1 ABC	ENTER	
# Enter the Titration End point Peak number for calculating the results			
Derivative Threshold (1-32000) :	5 MNO	0 -./	0 -./
Derivative Threshold (1-32000) : 500	ENTER		

Conc. (0.00001-99.99999) :

Conc. (0.00001-99.99999) : 0.100

1)BlankVolume 2)PrimaryDose 3)None :

Constant Code (1-20) :

Const 1 (0.0001-999999.99) :

Constant (0.0001-999999.99) : 1

Print Report 1)Full 2)Select 3)Condensed 4)None : 4

Titration Report 1)Y 2)N : 1

Titration Data Table 1)Y 2)N : 2

Titration mV/ μ l Graph 1)Y 2)N : 1

Titration Derivative Graph 1)Y 2)N : 2

Titration Time/ μ l Graph 1)Y 2)N : 2

Titration 2nd Deriv Graph 1)Y 2)N : 2

Method No. 2 Stored

After 4-5 seconds

1)Run 2)View :

• **TITRATION BY pH CUTOFF (END) VALUE.:-**

METHOD NO.	: 3	TITRANT NAME	: NAOH
SAMPLE NAME	: HCl	TITR. TYPE	: Acid-Base Aqueous
MEAS. TYPE	: Potentiometric	P.D.CRITERIA	: 0 mV
TITR. MODE	: Cutoff pH	P.D.DOSE	: 100 µl
P.D.VOLUME	: 500 µl	STIR TIME	: 5 sec
P.D.INTERVAL	: 5 sec	EQ. FACTOR	: 5
INITIAL DOSE	: 100 µl	CTRL. LIMIT	: 4.00 pH
VOLUME LIMIT	: 10000 µl	MIN.CTRL.DOSE	: 10 µl
END LIMIT	: 9.00 pH	BV/PD/None	: None
DERIV.THRESH.	: 500	CONSTANT	: 1
CONC. (mol/l)	: 0.100	T. DATA TABLE	: 2 (No)
CONSTANT CODE	: 1	DERIV. GRAPH	: 2 (No)
TITR. REPORT	: 1 (Yes)	2nd DERIV. GRAPH	: 2 (No)
mV/µl GRAPH	: 1 (Yes)		
Time/µl GRAPH	: 2 (No)		

1)Develop Method 2) Function 3)Config. : 7 4)Print 5) Display 6) Reprocess 7)Run	1 ABC
---	----------

1)Quick method 2)Full method :	2 DEF
--------------------------------	----------

1)Add/Edit method 2) Copy method 3) Skip method 4) View :	1 ABC
--	----------

Enter User Password:

Enter ADMIN or USER Password (ref. Page. 51)

Method No. (1-40) :	3 GHI	ENTER
---------------------	----------	-------

SAMPLE NAME :		ENTER
---------------	--	-------

TITRANT NAME :		ENTER
----------------	--	-------

Enter Sample and Titrant Name [refer Appendix-A for alphanumeric entries]

Measurement Type 1)Potentiometric 2)Voltametric :	1 ABC	ENTER
--	----------	-------

Titration Type 1)Aq. 2)Non-Aq. 3)Redox 4)Pr. 5)EDTA 6)Back 7)Photo 8)Thermo : 1	1 ABC	ENTER
--	----------	-------

Select titration mode 1)Incr. 2)Equil. 3)Cutoff(pH) 4)pHstat: 1	3 GHI	ENTER
--	----------	-------

Predispense Criteria (0-1000 mV) : 0	0 -./	ENTER
--------------------------------------	----------	-------

Predispense Volume (0-95000 µl) : 0

5	0	0
MNO	-./	-./

Predispense Volume (0-95000 µl) : 500

ENTER

Predispense Dose (5-1000) :

1	0	0
ABC	-./	-./

Predispense Dose (5-1000) : 100

ENTER

Predispense Interval (0-999 sec) : 0

5
MNO

ENTER

Stir Time (0-999 sec) : 5

5
MNO

ENTER

Initial Dose (5-2000 µl) : 100

1	0	0
ABC	-./	-./

Initial Dose (5-2000 µl) : 100

ENTER

Eq. Factor(1-99) : 5

5
MNO

ENTER

Volume Limit (Max. : 99999 µl) :

1	0	0
ABC	-./	-./

0	0
-./	-./

Volume Limit (Max. : 99999 µl) : 10000

ENTER

Control Limit (0-14 pH) :

4
JKL

.

0	0
-./	-./

Control Limit (0-14 pH) : 4.00

ENTER

End Limit (0-14 pH) :

9
YZ_

.

0	0
-./	-./

End Limit (0-14 pH) : 9.00

ENTER

Minimum Control Dose (5-500 µl) :

1	0
ABC	-./

Minimum Control Dose (5-500 µl) : 10

ENTER

Derivative Threshold (1-32000) :	5 MNO	0 -./	0 -./
Derivative Threshold (1-32000) : 500			ENTER
Conc. (0.00001-99.99999) :	0 -./	.	1 ABC
0 -./	0 -./	Conc. (0.00001-99.99999) : 0.100	ENTER
1)BlankVolume 2)PrimaryDose 3)None :	3 GHI		ENTER
Constant Code (1-20) :	1 ABC		ENTER
Const 1 (0.0001-999999.99) :			1 ABC
Const 1 (0.0001-999999.99) : 1			ENTER
Print Report 1)Full 2)Select 3)Condensed 4)None : 4	2 DEF		ENTER
Titration Report 1)Y 2)N : 1	1 ABC		ENTER
Titration Data Table 1)Y 2)N : 2	2 DEF		ENTER
Titration mV/ul Graph 1)Y 2)N : 1	1 ABC		ENTER
Titration Derivative Graph 1)Y 2)N : 2	2 DEF		ENTER
Titration Time/ul Graph 1)Y 2)N : 2	2 DEF		ENTER
Titration 2nd Deriv Graph 1)Y 2)N : 2	2 DEF		ENTER
Method No. 3 Stored			After 4-5 seconds
1)Run 2)View :			

• **TITRATION BY PH STAT METHOD: FOR ENZYME KINETIC STUDY.**

e.g. say pH to be maintained is 8 pH.±0.02pH

METHOD NO.	: 4	TITRANT NAME	: NAOH
SAMPLE NAME	: ENZYME	TITR. TYPE	: Acid-Base Aqueous
MEAS. TYPE	: Potentiometric	P.D.CRITERIA	: 0 mV
TITR. MODE	: pH stat	STIR TIME	: 5 sec
P.D.VOLUME	: 0 µl	STAT INTV.	: 5 sec
INITIAL DOSE	: 100 µl	CTRL. LIMIT	: 8.00 pH
VOLUME LIMIT	: 10000 µl	STAT RUN TIME	: 90 sec
END LIMIT	: 8.02 pH	VTIME 2	: 90 sec
VTIME 1	: 30 sec	MIN.CTRL DOSE	: 5 µl
FACTOR	: 1.000	BV/PD/None	: None
DERIV.THRESH.	: 50	CONSTANT	: 1
CONC. (mol/l)	: 0.100		
CONSTANT CODE	: 8		
PRINT REPORT	: Full		

Note : To understand the meaning of each parameter, please refer Appendix F-H.

1)Develop Method 2) Function 3)Config. : 7 4)Print 5) Display 6) Reprocess 7)Run	1
	ABC

1)Quick method 2)Full method :	2
	DEF

1)Add/Edit method 2) Copy method 3) Skip method 4) View :	1
	ABC

Enter User Password:	Enter ADMIN or USER Password (ref. Page. 51)
----------------------	--

Method No. (1-40) :	4	ENTER
	JKL	

SAMPLE NAME :	ENTER
---------------	-------

TITRANT NAME :	ENTER
----------------	-------

Enter Sample and Titrant Name [refer Appendix-A for alphanumeric entries]

Measurement Type 1)Potentiometric 2)Voltametric :	1	ENTER
	ABC	

Titration Type 1)Aq. 2)Non-Aq. 3)Redox 4)Pr. 5)EDTA 6)Back 7)Photo 8)Thermo : 1	1	ENTER
	ABC	

Select titration mode 1)Incr. 2)Equil. 3)Cutoff(pH) 4)pHstat: 4	4	ENTER
	JKL	

Predispense Criteria (0-1000 mV) : 0	0	ENTER
	- ./	

Predispense Volume (0-95000 µl) : 0 ENTER

Stir Time (0-999 sec) : 5 5 ENTER
MNO

Initial Dose (5-2000 µl) : 100 1 0 0
ABC -./ -./

Initial Dose (5-2000 µl) : 100 ENTER

Stat Interval (1-99 sec) : 5 5 ENTER
MNO

Volume Limit (Max. : 99999 µl) : 1 0 0
ABC -./ -./

0 0 Volume Limit (Max. : 99999 µl) : 10000 ENTER
-./ -./

Control Limit (0-14 pH) : 8 .
VWX

0 0 Control Limit (0-14 pH) : 8.00 ENTER
-./ -./

End Limit (0-14 pH) : 8 .
VWX

0 2 End Limit (0-14 pH) : 8.02 ENTER
-./ DEF

Stat Run Time (1-9999 sec) : 9 0
YZ_ -./

Stat Run Time (1-9999 sec) : 90 ENTER

VTime 1 (1-9999 sec) : 6 0
PQR -./

VTime 2 (1-9999 sec) : 9 0
YZ_ -./

Factor (0.001-10000) : 1 . 0
ABC -./

0 0 Factor (0.001-10000) : 1.000 ENTER
-./ -./

Minimum Control Dose (5-500 µl) : 1 0
ABC -./

Minimum Control Dose (5-500 μ l) : 10 ENTER

Derivative Threshold (1-32000) : 2 0 0
DEF -./ -./

Derivative Threshold (1-32000) : 200 ENTER

Conc. (0.00001-99.99999) : 0 . 1
-./ -./ ABC

0 0 Conc. (0.00001-99.99999) : 0.100 ENTER
-./ -./

1)BlankVolume 2)PrimaryDose 3)None : 3 ENTER
GHI

Constant Code (1-20) : 8 ENTER
VWX

Print Report 1)Full 2)Select 3)Condensed
4)None : 4 ENTER
JKL

Method No. 4 Stored

After 4-5 seconds

1)Run 2)View :

• **TITRATION BY MULTIPLE END POINT METHOD:**

METHOD NO.	: 5	TITRANT NAME	: AgNO ₃
SAMPLE NAME	: HALIDES	TITR. TYPE	: Precipitation
MEAS. TYPE	: Potentiometric	P.D.CRITERIA	: 0 mV
TITR. MODE	: Incremental	STIR TIME	: 5 sec
P.D.VOLUME	: 0 µl	INTERVAL	: 5 sec
INITIAL DOSE	: 100 µl	TITR. END CRIT.	: End Point
VOLUME LIMIT	: 5000 µl	CALCULATION BY	: ALL
NO. OF EP	: 3	BV/PD/None	: None
DERIV.THRESH.	: 20	CONSTANT 1	: 16.9
CONC. (mol/l)	: 0.100	CONSTANT 2	: 11.9
CONSTANT CODE	: 13	CONSTANT 3	: 3.55
PRINT REPORT	: None		

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run	1
	ABC

1)Quick method 2)Full method :	2
	DEF

1)Add/Edit method 2) Copy method 3) Skip method 4) View :	1
	ABC

Enter User Password: Enter ADMIN or USER Password (ref. Page. 51)

Method No. (1-40) :	5	ENTER
	MNO	

If method is already present , Display reads :

Method exists !! Edit 1)Y 2)N :	1
	ABC

SAMPLE NAME :	ENTER
---------------	-------

TITRANT NAME :	ENTER
----------------	-------

Enter Sample and Titrant Name [refer Appendix-A for alphanumeric entries]

Measurement Type 1)Potentiometric 2)Voltametric :	1	ENTER
	ABC	

Titration Type 1)Aq. 2)Non-Aq. 3)Redox 4)Pr. 5)EDTA 6)Back 7)Photo 8)Thermo : 1	4	ENTER
	JKL	

Select titration mode 1)Incr. 2)Equil. 3)Cutoff(pH) 4)pHstat: 2	1	ENTER
	ABC	

Predispense Criteria (0-1000 mV) : 0
- ./ ENTER

Predispense Volume (0-95000 µl) : 0
- ./ ENTER

Stir Time (0-999 sec) : 5
MNO ENTER

Initial Dose (5-2000 µl) : 1 0 0
ABC - ./ - ./

Initial Dose (5-2000 µl) : 100 ENTER

Interval (1-999 sec) : 5
MNO ENTER

Volume Limit (Max. : 99999 µl) : 5000 5 0
MNO - ./

0 0 Volume Limit (Max. : 99999 µl) : 5000 ENTER
- ./ - ./

Titration End Criteria
1] EndPoint (EP) 2] VolumeLimit (VL) : 1 1
ABC ENTER

Number of EP (1-9) : 3
GHI ENTER

Calculation by
1)1st 2)Largest 3>Last 4)All 5)Select : 1 4
JKL ENTER

Auto Eval 1)Y 2) N : 1 ENTER

Derivative Threshold (1-32000) : 2 0
DEF - ./

Derivative Threshold (1-32000) : 20 ENTER

Conc. (0.00001-99.99999) : 0 . 1
- ./ ABC

0 0 Conc. (0.00001-99.99999) : 0.100 ENTER
- ./ - ./

1)BlankVolume 2)PrimaryDose 3)None : 3
GHI ENTER

1 3
ABC GHI

Constant Code (1-20) :

ENTER

Const 1 (0.0001-999999.9) :

Say – 16.9, Press keys 1,6, dot & 9

Const 1 (0.0001-999999.9) : 16.9

ENTER

Const 2 (0.0001-999999.9) :

Say – 11.9, Press keys 1,1, dot & 9

Const 2 (0.0001-999999.9) : 11.9

ENTER

Const 3 (0.0001-999999.9) :

Say – 3.55, Press keys 3, dot, 5 & 5

Const 3 (0.0001-999999.9) : 3.55

ENTER

Print Report 1)Full 2)Select 3)Condensed 4)None : 4

4
JKL

ENTER

Method No. 5 Stored

After 4-5 seconds

1)Run 2)View :

TITRATION METHOD FOR BLANK & BACK TITRATION:-

Back-titrations may be used in volumetric analysis where direct titrations cannot be used for technical reasons such as Sample solubility, Trace of impurities, and slow reaction. For example, when a compound is not soluble in water, or where it contains impurities, which interfere with direct titration, an excess of the titrant may be added, and the excess then determined by titration (i.e. by back titration). The volume required for sample is determined by the difference between the two titrants endpoint volume.

Firstly Perform determine the Blank Volume for Excess Titrant say 20mL of 0.1M NaOH by titration with 0.1M HCL. The volume required for endpoint (EP volume) is entered as Blank Volume in the Back titration method during the actual BACK titration

METHOD NO.	: 6	(To estimate Blank Volume)	
SAMPLE NAME	: NaOH	TITRANT NAME	: HCl
MEAS. TYPE	: Potentiometric	TITR. TYPE	: Acid-Base Aqueous
TITR. MODE	: Equilibrium	P.D.CRITERIA	: 0 mV
P.D.VOLUME	: 18000 µl	P.D.DOSE	: 500 µl
P.D.INTERVAL	: 5 sec	STIR TIME	: 30 sec
INITIAL DOSE	: 200 µl	EQ. FACTOR	: 6 sec
VOLUME LIMIT	: 22000 µl	TITR. END CRIT.	: End Point
NO. OF EP	: 1	CALCULATION BY	: 1 st
MIN.CONT.DOSE	: 10	DERIV.THRESH.	: 300
		CONC. (mol/l)	: 0.100
BV/PD/None	: None	CONSTANT CODE	: 8
PRINT REPORT	: None		

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

1
ABC

1)Quick method 2)Full method :

2
DEF

1)Add/Edit method 2) Copy method
3) Skip method 4) View :

1
ABC

Enter User Password:

Enter ADMIN or USER
Password (ref. Page. 51)

Method No. (1-40) :

6
PQR

ENTER

If method is already present , Display reads :

Method exists !! Edit 1)Y 2)N :

1
ABC

SAMPLE NAME :

ENTER

TITRANT NAME :

ENTER

Enter Sample and Titrant Name [refer Appendix-A for alphanumeric entries]

Measurement Type
1)Potentiometric 2)Voltametric :

1
ABC

ENTER

Titration Type 1)Aq. 2)Non-Aq. 3)Redox
4)Pr. 5)EDTA 6)Back 7)Photo 8)Thermo : 1

1
ABC

ENTER

Select titration mode
1)Incr. 2)Equil. 3)Cutoff(pH) 4)pHstat: 2

2	ENTER
DEF	

Predispense Criteria (0-1000 mV) : 0

0	ENTER
-./	

Predispense Volume (0-95000 µl) : 0

Enter 1,8,0,0,0 and

ENTER

Predispense Volume (0-95000 µl) : 18000

ENTER

Predispense Dose (5-1000) :

5	0	0
MNO	-./	-./

Predispense Dose (5-1000) : 500

ENTER

Predispense Interval (0-999 sec) : 0

5	ENTER
MNO	

Stir Time (0-999 sec) :

3	0
JKL	-./

Stir Time (0-999 sec) : 30

ENTER

Initial Dose (5-2000 µl) :

2	0	0
DEF	-./	-./

Initial Dose (5-2000 µl) : 200

ENTER

Eq. Factor(1-99) :

6	ENTER
PQR	

Volume Limit (Max. : 99999 µl) :

2	2	0
DEF	DEF	-./

0	0	Volume Limit (Max. : 99999 µl) : 22000	ENTER
-./	-./		

Titration End Criteria
1] EndPoint (EP) 2] VolumeLimit (VL) : 1

1	ENTER
ABC	

Number of EP (1-9) :

1	ENTER
ABC	

Calculation by
 1)1st 2)Largest 3>Last 4)All 5)Select : 1

1
 ABC ENTER

Minimum Control Dose (5-500 µl) :

1 0
 ABC -./

Minimum Control Dose (5-500 µl) : 10

ENTER

Derivative Threshold (1-32000) :

3 0 0
 GHI -./ -./

Derivative Threshold (1-32000) : 300

ENTER

Conc. (0.00001-99.99999) :

0 . 1
 -./ ABC

0 0 Conc. (0.0001-99.99999) : 0.100
 -./ -./ ENTER

1)BlankVolume 2)PrimaryDose 3)None :

3
 GHI ENTER

Constant Code (1-20) :

8
 VWX ENTER

Constant : 1

Print Report 1)Full 2)Select 3)Condensed
 4)None : 4

4
 JKL ENTER

Method No. 6 Stored

After 4-5 seconds

1)Run 2)View :

Run this Method to obtain the volume required for the titrant to be added as excess in the Back Titration. (Refer to RUN part in this manual)

• **METHOD FOR BACK TITRATION:-**

METHOD NO.	: 7 (BACK Titration)	TITRANT NAME	: HCl
SAMPLE NAME	: Sample + NaOH	TITR. TYPE	: Back
MEAS. TYPE	: Potentiometric	P.D.CRITERIA	: 0 mV
TITR. MODE	: Equilibrium	P.D.DOSE	: 500 µl
P.D.VOLUME	: 10000	STIR TIME	: 30 sec
P.D.INTERVAL	: 5 sec	EQ. FACTOR	: 6 sec
INITIAL DOSE	: 200 µl	TITR. END CRIT.	: End Point
VOLUME LIMIT	: 22,000	CALCULATION BY	: 1 st
NO. OF EP	: 1	DERIV.THRESH.	: 300
MIN.CONT.DOSE	: 10	CONC. (mol/l)	: 0.100
BV/PD/None	: 19.5	CONSTANT CODE	: 13
PRINT REPORT	: None		

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run	1
	ABC

1)Quick method 2)Full method :	2
	DEF

1)Add/Edit method 2) Copy method 3) Skip method 4) View :	1
	ABC

Enter User Password: Enter ADMIN or USER Password (ref. Page. 51)

Method No. (1-40) :	7	ENTER
	STU	

If method is already present, Display reads :

Method exists !! Edit 1)Y 2)N :	1
	ABC

SAMPLE NAME :	ENTER
---------------	-------

TITRANT NAME :	ENTER
----------------	-------

Enter Sample and Titrant Name [refer Appendix-A for alphanumeric entries]

Measurement Type 1)Potentiometric 2)Voltametric :	1	ENTER
	ABC	

Titration Type 1)Aq. 2)Non-Aq. 3)Redox 4)Pr. 5)EDTA 6)Back 7)Photo 8)Thermo : 1	6	ENTER
	PQR	

Select titration mode 1)Incr. 2)Equil. 3)Cutoff(pH) 4)pHstat: 2	2	ENTER
	DEF	

Predispense Criteria (0-1000 mV) : 0	0	ENTER
	-./	

Predispense Volume (0-95000 µl) : 0 Enter 10000 and ENTER

Predispense Volume (0-95000 µl) : 10000 ENTER

Predispense Dose (5-1000) : 5 0 0
MNO - ./ - ./

Predispense Dose (5-1000) : 500 ENTER

Predispense Interval (0-999 sec) : 0 5 ENTER
MNO

Stir Time (0-999 sec) : 3 0
JKL - ./

Stir Time (0-999 sec) : 30 ENTER

Initial Dose (5-2000 µl) : 2 0 0
DEF - ./ - ./

Initial Dose (5-2000 µl) : 200 ENTER

Eq. Factor(1-99) : 6 ENTER
PQR

Here the Volume limit is entered

Volume Limit (Max. : 99999 µl) : 2 2 0
DEF DEF - ./

0 0 Volume Limit (Max. : 99999 µl) : 22000 ENTER
- ./ - ./

Titration End Criteria
1] EndPoint (EP) 2] VolumeLimit (VL) : 1 1 ENTER
ABC

Number of EP (1-9) : 1 ENTER
ABC

Calculation by
1)1st 2)Largest 3>Last 4)All 5)Select : 1 1 ENTER
ABC

Minimum Control Dose (5-500 µl) : 1 0
ABC - ./

Minimum Control Dose (5-500 µl) : 10 ENTER

Derivative Threshold (1-32000) : 3 0 0
GHI -./ -./

Derivative Threshold (1-32000) : 300 ENTER

Conc. (0.00001-99.99999) : 0 . 1
-./ -./ ABC

0 0 Conc. (0.0001-99.99999) : 0.100 ENTER
-./ -./

1)BlankVolume 2)PrimaryDose 3)None : 1 ENTER
ABC

Blank Volume (0.0001-99.9999ml) :

Enter the Blank volume determined for the Excess NAOH added from the 1st titration

Say 19.525ml

Volume Factor (0.100-5.000) : 1.0000

Constant Code (1-20) : 1 3 ENTER
ABC GHI

Const. 1 (0.0001-999999.9) : 1.000

Now select the appropriate Result calculation Constant and Enter the Value (refer the Previous Method generation examples)

Print Report 1)Full 2)Select 3)Condensed 4)None : 4 ENTER
JKL

Method No. 7 Stored

After 4-5 seconds

1)Run 2)View :

Entering Method for TAN & TBN (ASTM D664, ASTM D2896 & ASTM D4739)

Note : Select following Titration methods numbers for TAN & TBN

- **Method No. 39** for **ASTM D2896**, which is indicated on printout.
- **Method No. 38 & 40** are respectively for **ASTM D664** and **D4739** for which buffer values should be entered as following, the same is indicated on printout.

Buffer mV Value (1-3200) :	2	0
	DEF	- ./

0	Buffer mV Value (1-3200) : 200	ENTER
- ./		

Sign 1)+ 2)- :	1	ENTER
	ABC	

Titration Trend 1)+ 2)- 3)None :	1	ENTER
	ABC	

- **Viewing the Titration Methods :** All the 50 methods can be scanned for checking or modifying the parameters individually.

Say Method No.1 has to be viewed before performing analysis.

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run	1 ABC	
1)Quick method 2)Full method :	2 DEF	
1)Add/Edit method 2) Copy method 3) Skip method 4) View :	4 JKL	ENTER
Method No. (1-50) :	1 ABC	ENTER
SAMPLE NAME : HCl		ENTER
TITRANT NAME : NAOH		ENTER
Measurement Type : Potentiometric		ENTER
Type : Acid-Base Aqueous Mode : Incremental		ENTER
Predispense Criteria (mV) : 0 Predispense Volume (μl) : 500		ENTER
Predispense Dose (μl) : 100 Predispense Interval (sec) : 5		ENTER
Stir Time (sec) : 5 Initial Dose (μl) : 100		ENTER
Interval (sec) : 5 VolLimit (μl) : 10000		ENTER
Titr. End Crit. : EndPoint No. of EP : 1 Cal. By : 1st		ENTER
Deriv. Threshold : 50		ENTER
Concentration : 0.10000		ENTER
BV/PD/None : None		ENTER
Constant Code : 1 Constant : 63.0000		ENTER

Print Report : None


ENTER

**1)Add/Edit method 2) Copy method
3) Skip method 4) View :**

NOTE : 50 titration methods can be scanned independently. If the method scanned is not present or earlier entered, display shows :

Method does not exist !!


- **Copying the Titration Method :** Say Method No.1 has to be copied from Method No. 41. Follow the below mentioned instructions




• **Skipping the Titration Method :**

Say Method No.1 has to be skipped so memory can be made free to add new method.

a) To skip individual method :



b) To skip all methods (except standard methods 41-50) :



FUNCTION PARAMETERS

1) pH Calibration :- for pH measurement or pH titrations.

Standardization is done at 7 pH and Calibration can be done for acid and base buffer. Calibration can be performed for 3 point with 2 buffer (one on acidic side and other on basic side). If the calibration is not to be done, Press “**Escape**” key to quit.

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run	2
	DEF

a) Calibration of pH electrode :- Say 7 pH , 4 pH and 9pH buffer

Remove the pH electrode from its storage buffer and rinse it with distilled water. Clean the electrode soak dry with tissue paper and put the electrode in 7 pH buffer.

1)pH Calibration 2)mV scan 3)Prime 4)Burette Change 5)Manual Titra :	1	ENTER
	ABC	

Enter User Password:

Enter ADMIN or USER
Password (ref. Page. 51)

1)Calibrate 2)Clear Calibration	1	ENTER
	ABC	

Dip the electrode in 7.00 pH buffer for standardization,

Enter if stable mV (std) mV = 0.00	ENTER
---------------------------------------	-------

System checks the Iso-potential value in ± 30 mV limit and display reads,

Offset : 1.2 mV Standardisation Over ...

After 4-5 seconds

Note : 1) The offset shown here is now memorized for future measurements/analysis.
2) If the offset is out of the range (more than +/- 30 mV), the display shows

Incorrect Buffer or Faulty Electrode

Check the buffer or change the electrode and proceed further.

To complete the calibration process use buffer of acid side first and then basic side

Remove the pH electrode from 7 pH buffer and rinse it with distilled water. Sock dry the electrode with tissue paper to remove the excess DI water and put the electrode in Acidic buffer (say 4 pH).

Enter pH Value (0-6.5 pH) :	4	.
	JKL	

0	0	Buffer 1 (0-6.5 pH) : 4.00	ENTER
-./	-./		

Enter if stable mV (std) mV = +177.0	T : 25.0°	ENTER
---	-----------	-------

Remove the pH electrode from 4 pH buffer and rinse it thoroughly with distilled water. Sock dry the electrode with tissue paper to remove the excess DI water and put the electrode in Base buffer (say 9 pH).

Buffer 2 (7.5-14 pH) :	1 ABC	0 -./	•
------------------------	----------	----------	---

0 -./	0 -./	Buffer 2 (7.5-10 pH) : 10.00	ENTER
----------	----------	------------------------------	-------

Enter if stable mV (std) mV = +177.0	T : 25.0°	ENTER
---	-----------	-------

Calibration Over ...

After 4-5 seconds

Data Transfer to PC

After 4-5 seconds screen returns to main menu.

The electrode parameters (offset, slope1, slope2) are calculated, automatically stored and printed in 'Calibration report' along with Calibration date and time (if printer is connected). The slope limit provided in the instrument is 80% to 120%)

c) Resetting the Calibration data to Default Values

If instrument calibration done by buffers is not required and needs to be resetted press '2' , the system takes default values with slope (100%) and offset correction (0.00mV).

1)Calibrate 2)Clear Calibration	2 DEF	ENTER
---------------------------------	----------	-------

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run


After 4-5 seconds screen returns to main menu.

‘ Now the instrument is ready for pH measurement use.’

2) Using instrument as pH / mV indicator (mV scan)

For pH measurement ensure that the instrument is properly calibrated before measurement.


Place the beaker containing sample whose pH is to be measured. Insert the clean and dried electrode such that it gets immersed in the sample. Now start the stirrer and proceed as following



3) Priming To rinse the complete liquid path by titrant and to remove air/air bubbles trapped from tubings .

In this process burette will aspire entered volume from titrant bottle and dispense it in reaction vessel or beaker. In case air bubbles in burette are noticed, tighten teflon tubing connectors. Initially PRIME function can be repeated 2 to 3 times to clean system.

The PRIME is also used for adding fixed volume as per program parameter = Primary Dose (**P.Dose**) is added.



Wait

After 4-5 seconds & Priming


Priming Over.

Screen returns to main menu.

NOTE : Max.volume that can be primed is 65 ml. Min.volume is 1/1000th burette.

The Priming Volume should always be in multiple of 1/1000th of the burette volume.

4) Replacing the burette assembly



Note: Burette Factor is for volume correction related to dispensing. The factor is multiplied to the End point/equivalence point Volume obtained for correcting the results.

The instruments automatically restarts and confirms the detection and returns to Main Menu.

5) Performing Manual Titration

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

2
DEF

1)pH Calibration 2)mV scan 3)Prime
4)Burette Change 5)Manual Titra:

5
MNO

Titration Dose (5-10000 ul) :

1	0	0
ABC	-./	-./

Titration Dose (5-10000 ul) : 100

ENTER

Wait

After refilling burette

Press ENTER to Continue / ESC to Exit

Press ESC to Exit Or
After pressing ENTER

Disp.Vol. :	25.0	BV 10	TITRA
pH : xx.yy	mV : aaa.bb	00 : 00	: 01


Press ENTER to dose or
ESC to exit

pH, mV values and Volume dispensed are updated after every dose is added by pressing the "ENTER" key.

After "ESC" key is pressed, system initializes and data table report from main menu can be printed.


CONFIG PARAMETERS

1) Setting RTC (Clock)




After entering the values (new time, date are stored), screen returns to main menu.

2) Balance Connectivity : for automatic sample weight transfer to instrument.




a) For Manual Sample weight entry using Keys



Screen returns to main menu.

b) For Sample weight entry direct from Balance (Weighing Scale) RS232




Screen returns to main menu.

Balance Setting :-


For Sartorius Balance : Set the Comm Port Properties, Baud – 1200, Data bits – 8, Parity – Odd, Stop bit – 1.

For functional details refer the titration Run Display


3) Temperature Compensation



a) For Manual Temperature compensation




Enter the sample temperature after measuring it with external thermometer (25)





Screen returns to main menu.

b) For Automatic Temperature compensation



Screen returns to main menu.




4) Printer Selection




To select LX800



5) Password Entry



Min. 4 digit Max. 8 digit



Admin Password saved

To enter User Password

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

3

GHI

1)Clock/Date 2)Balance 3)Temp. 4)Printer
5)PW 6)Data Xfer 7) Resol. 8) Tit. Temp. Fct.

5

MNO

1)Admin P/W 2)User P/W :

To change

2

DEF

Enter Admin Password :

Set User Password :

Min. 4 digit Max. 8 digit

Repeat Password :

User Password saved

6) DATA TRANSFER

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

3

GHI

1)Clock/Date 2)Balance 3)Temp. 4)Printer
5)PW 6)Data Xfer 7) Resol. 8) Tit. Temp. Fct.

6

PQR

1) Printer 2) PC

1

ABC

For Data Transfer to Printer

OR

2

DEF

For Data Transfer to PC
(HyperTerminal)

After the data is printed / transferred to PC, Display returns to

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

7) Resolution selection

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run	3 GHI
---	----------

1)Clock/Date 2)Balance 3)Temp. 4)Printer 5)PW 6)Data Xfer 7) Resol. 8) Tit. Temp. Fct.	7 STU
---	----------

Result Resolution : 1)1 2)2 3)3 4)4

Enter the required result resolution, Say 4

Press

4 JKL

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run

8) Titrant Temp Factor

1)Develop Method 2) Function 3)Config. 4)Print 5) Display 6) Reprocess 7)Run	3 GHI
---	----------

1)Clock/Date 2)Balance 3)Temp. 4)Printer 5)PW 6)Data Xfer 7) Resol. 8) Tit. Temp. Fct.	6 PQR
---	----------

Titrant Temp Factor (0.900-1.100):	1 ABC	ENTER
------------------------------------	----------	-------

This factor will get considered in calculations.

PRINT

1) Printing Report

After the titration is over, the result is stored in non-volatile memory, which can be printed for report generation. The Hardcopy of titration can be taken for : i) Result
ii) Data iii) Graph (mV v/s μl or 1st derivative of mV v/s μl) iv) Condensed Report

Care has been taken for printing in DRAFT mode. Result gives method and mode used, titration number, sample name, titrant name, run identification number and formula constants. It gives result report as end point ml with desired result units depending on the constant selection. Data gives stepwise μl reading with corresponding mV, S.No. and change in mV with each step of dosing. Graphical representation if data table is graph. User can call any one of these reports as per his requirement. Data can be transferred and stored in PC for future reference.

If the second element is 0 (i.e. this condition occurs only after resetting the system), graph cannot be printed, error message is displayed.

ERROR 2

If number of doses exceeds 199, during any titration run, graph cannot be printed, error message is displayed.

ERROR 2

If first derivative value is out of range (-32000 to 32000), during any titration run, graph cannot be printed, error message is displayed.

ERROR 3

NOTE : Auto Incremented Run number is printed on the right side of every printout. If the volume reqd. is negative (blank/back volume exceeds titration volume) the result is printed and displayed **NEGATIVE**.

• Printing Last Titration Result

Ensure that the calculation parameters are entered before taking hardcopy. If any of calculation parameter is altered after 1st titration run report, '(•)' is printed near Run No.

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

4

JKL

The Print function performs its operation as per configuration done in Data Xfer. For Data Transfer mode selected to 'Printer' display shows,

1)Result 2)Data 3)Condensed 4) Statistics
5) Cal. Repo 6) Method 7) Graph

Press

1

ABC

To take Result Printout


Titration Report is printed and screen returns to main menu.

For Data Transfer mode selected to 'PC' display shows,

Downloading to PC


Titration Data is transferred to PC and screen returns to main menu.

- **Printing Data Table**




Titration Data Table is printed and screen returns to main menu.

- **Printing the Titration Analysis Condensed Report** Follow the below mentioned instructions :



Titration Analysis Condensed Report is printed and screen returns to main menu.

- **Printing the Titration Statistics Report :** Follow the below mentioned instructions to print statistics report for 3 readings without skipping any of them.




Titration Statistics Report is printed and screen returns to main menu.

For No data it will display


No data for statistics

- **Printing pH Calibration Report :** Follow below mentioned instructions :



pH Calibration Report with Buffer values and Slope is printed along with the last calibration date is printed and screen returns to main menu.



• **Printing Method Parameters** Follow below mentioned instructions





Parameter Report is printed along with empty methods & screen returns to main menu.





Select Prog. No. Say 1



If printer is not connected / Ready



• **Printing the mV v/s μ l Graph Report**



Titration Graph Report is printed and screen returns to main menu

Similarly other graphs can be print.

Note : Printouts of Graphs are not available if Data Transfer mode selected to PC.

DISPLAY

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

5
MNO

1) Displaying Titration Result Follow below mentioned instructions

1) Result 2) Statistics 3) Cal. Data

1
ABC

Vol. Reqd : 0.9425 ml
Result : 0.1924 Molarity (T)

After 4-5 sec screen

Result - 1)Print 2)Reprocess 3)Display
4)Display Stat. 5)Run Last Method :

2) Displaying the Titration Statistics Report : Follow the below mentioned instructions to display statistics report for 3 readings without skipping any of them.

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

5
MNO

1) Result 2) Statistics 3) Cal. Data

2
DEF

Readings (2-10) :

3
GHI

ENTER

Skip Reading (0-10) :

0
./

ENTER

Titration Statistics are displayed and screen returns to main menu.

3) Displaying the Cal. Data :

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

5
MNO

1) Result 2) Statistics 3) Cal. Data

3
DEF

Offset (mV): X DD/MM/YY
Slope 1 : ABC% Slope 2 : DEF%

REPROCESSING THE RESULT FOR DIFFERENT RESULT UNITS :

Derivative Threshold (1-32000) :	2 DEF	0 -./
----------------------------------	----------	----------

Derivative Threshold (1-32000) : 20	ENTER
-------------------------------------	-------

Calculation by 1)1 st 2)Largest 3>Last 4)All 5)Select : 1	1 ABC	ENTER
---	----------	-------

Conc. (0.00001-99.99999) :	0 -./	.	1 ABC
----------------------------	----------	---	----------

0 -./	0 -./	Conc. (0.00001-99.99999) : 0.100	ENTER
----------	----------	----------------------------------	-------

1)BlankVolume 2)PrimaryDose 3)None :	3 GHI	ENTER
--------------------------------------	----------	-------

Constant Code (1-8) :	1 ABC	ENTER
-----------------------	----------	-------

Constant (0.0001-999999.99) :	6 PQR	3 GHI
-------------------------------	----------	----------


Constant (0.0001-999999.99) : 63	ENTER
----------------------------------	-------

Sample Vol. (0.0001 – 999.99999ml) :


Vol. Reqd : 0.9425 ml Result : 0.1924	Molarity (s)
--	--------------

After 4-5 sec Screen returns to Main Menu


RUNNING METHOD : Follow below mentioned instructions :





a) Run last method




If last method is quick method, Press ENTER to acknowledge :



If last method is full method, Press ENTER to acknowledge :




Enter Identification No. [refer Appendix-A for alphanumeric entries]




NOTE :


If constant (1-20) for the method is related to weight and balance mode is MANUAL, screen shows :




If constant (1-20) for the method is related to weight and balance mode is AUTO, screen shows :



After pressing Print key of balance, screen changes to :



After pressing ENTER key to acknowledge on instrument screen changes to :



Wait
00 : 00 : 18

Method No. 1
starts with the
mentioned para.

Disp.Vol. : XXX ul # RN BV 10 TITRA
25.0° mV : 00 : 00 :

The instrument displays volume and corresponding mV after every dosing. After detecting the Titration end criteria, the exact volume required for endpoint and result is displayed :

Wait

Vol. Reqd : 0.9425 ml
Result : 1.0025 Molarity (T)

After displaying the result for 4-5 seconds, the selected report (Titration Analysis or Data Table or Graph Report) printout in method are printed and the screen changes to :

Result - 1)Print 2)Reprocess 3)Display
4)Display Stat. 5)Run Last Method :

b) Run new method

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

7
STU

1)Last Method 2)New Method :

2
DEF

Method No. (1-50) :

2
DEF

ENTER

Press ENTER to acknowledge and continue the run or ESC to exit to main menu.

IDENTIFICATION NO :

ENTER

Enter Identification No. [refer Appendix-A for alphanumeric entries]

Sample Volume (0.0001-99.9999 ml) :

1
ABC

ENTER

Running Method No. 2

Wait
00 : 00 : 24

Method No. 2 starts with the mentioned para.

Disp.Vol. : XXX ul # RN BV 10 TITRA
25.0° mV : 00 : 00 : 38

The instrument displays volume and corresponding mV after every dosing. After detecting the Titration end criteria, the exact volume required for endpoint and result is displayed :

Vol. Reqd : 0.9425 ml
Result : 1.0025 Molarity (T)

After displaying the result for 4-5 seconds, the selected report (Titration Analysis or Data Table or Graph Report) printout in method are printed and the screen changes to :

Result - 1)Print 2)Reprocess 3)Display
4)Display Stat. 5)Run Last Method :

c) Run last quick method

1)Develop Method 2) Function 3)Config.
4)Print 5) Display 6) Reprocess 7)Run

7
STU

1)Last Method 2)New Method :

1
ABC

If last method is quick method, Press ENTER to acknowledge, ESC to exit :

Run Quick Method
Press ENTER to continue / ESC to exit.

ENTER

1)Sample Weight 2)Sample Volume :

2
DEF

Sample Volume (0.0001-99.9999 ml) :

1
ABC

ENTER

Running Quick Method ...

Wait
00 : 00 : 18

Quick Method starts with the mentioned para.

Disp.Vol. : XXX ul # RN BV 10 TITRA +
25.0° mV : 00 : 00 :

Vol. Reqd : 0.9425 ml
1)Repeat Run 2)Result 3)Print :

- a) To repeat the same quick method with same parameters, select '1' option.
- b) To print the data table report for quick method, select '3' option.
- c) To get result for the quick method enter sample analysis parameters as following by selecting '2' option.

Calculation By:
1) 1st 2) Largest 3) Last 4) All 5) Select : 1

Conc. (0.00001-99.99999) :

0
-. /

.

1
ABC

0
-. /

0
-. /

Conc. (0.00001-99.99999) : 0.100

ENTER

1)BlankVolume 2)PrimaryDose 3)None :

3
GHI

ENTER

Constant Code (1-8) :

1
ABC

ENTER

Constant (0.0001-999999.99) : 1

1
ABC

ENTER

Sample Vol. (0.0001 – 999.99999 ml) :

1
ABC

Vol.Reqd : 0.9425ml R : 0.1924 MOLARITY(S)
1)Repeat Run 2)Reprocess 3)Print :

If the number of doses exceeds 200 or titration run volume exceeds 99.99 ml or after 2 doses Ctrl. Lt. for pH titration exceeds sample pH, error message is displayed, the titration run is aborted and suitably printed.

ERROR 1

SECTION 4

- **Result Calculation**
- **Entering Blank Volume / Primary Dose**
- **Entering Constant Related to Sample**
- **Automatic Determining and Storing of Molarity**
- **Maintenance and Cleaning of Burette**
- **Maintenance and Cleaning of Electrode**
- **Appendix A — K**

- **Result Calculation**

a) Direct Titration :-

The result is calculated according to the following formula :

$$\text{Result} = M1 * EP1 * C1 / S1$$

where,

EP1 = Titrant consumption up to the equivalence point in ml.

C1 = Concentration of titrant in mol./lit .

M1 = Calculation constant to convert the raw result into the desired unit. Refer Pg. No. 70

S1 = Weight (or Volume) of the sample in ml / gm.

b) Direct Titration with Blank or Primary dose:-

The result is calculated according to the following formula :

$$\text{Result} = M1 * (EP1+E0) * C1 / S1$$

where,

EP1 = Titrant consumption up to the equivalence point in ml.

E0 = Blank volume (subtracted) OR Primary dose (added) in ml.

C1 = Concentration of titrant in mol./lit.

M1 = Calculation constant to convert the raw result into the desired unit. It can be determined in accordance with Pg. No. 69

S1 = Weight (or Volume) of the sample in ml / gm.

• **Entering Blank Volume / Primary Dose**

a) **Blank Volume** : Blank Volume is the titrant volume required for neutralizing the known amount of solvent/distilled water taken for dissolving the sample. In certain titrations a small part of titrant is consumed by the solvent before actual sample endpoint is obtained. In such cases exact result requires correction by blank volume.

The blank volume gets subtracted from titrated volume.

e.g. say the 50 ml. Solvent / distilled water got neutralized by 100 µl (0.1 ml) of titrant.

1)BlankVolume 2)PrimaryDose 3)None :	1	ENTER
	ABC	

Blank Volume (0.0001-99.9999 ml) :	0	0	1
	-./	-./	ABC

0	0	Blank Volume (0.0001-99.9999 ml) : 0.100	ENTER
-./	-./		

Constant Code (1-20) :

b) **Primary Dose (PDose Volume)** : Primary Dose is the fixed Titrant volume added in the titration beaker when i) sample volume/weight taken in large size which will require end point volume more than the burette volume and ii) to reduce the titration time.

e.g. if the titrant volume required for neutralization is about 7 ml and the installed burette assembly is of 5ml capacity, then 5ml titrant can be directly dispensed using PRIME function (refer Page 50) and the same value should be entered in P.Dose and run the titration method to determine the end-point.

1)BlankVolume 2)PrimaryDose 3)None :	2	ENTER
	DEF	

PDose Volume (0.0001-99.9999 ml) :	5	.	0
	MNO		-./

0	0	PDose Volume (0.0001-99.9999 ml) : 5.000	ENTER
-./	-./		

Constant Code (1-20) :

NOTE : Last entered Blank Volume and Primary Dose are taken for result evaluation. Range for entering Blank Volume & Primary Dose is from 0.0001 to 99.9999 ml.

• **Entering CONSTANT related to Sample**



Enter Constant value by selecting appropriate code no [refer Appendix-B(Pg.75), J(Pg.86)] for automatic calculation of results in units (% , factor, molarity, ppm etc).

Code no 1 - 7 : are reserved and used for Liquid samples. Sample taken by volume (ml).

Code no 8 & 20 : has Constant value '1' and is to be used for End point titration when only titration end point volume is to be determined for liquid & solid samples respectively .






Code no 12-19 : are reserved and used for solid samples. Sample is taken by weight (gm).

e.g. Entering Constant value - 1 in Code no. 1 for Sample Volume of 1 ml. Constant code number is displayed on the bottom corner of the screen.



OR

e.g. Entering constant value - 0.2042 in Code no. 12 for 50.15 mg of Potassium Hydrogen Pthalate (KHP).



CONSTANT CODES

Suitable selection of constant code no. allows printout of the associated unit with result.

1) Constants (1-8) for sample taken by Volume (ml.)

CODE NO.	UNITS	FORMULA
1.	MOLARITY(S)	$Z1 / Z2$
2.	% VOL.	$(M*Z1) / (Z2*10 * \delta)$
3.	MG / L	$(M*1000 *Z1) / Z2$
4.	PPM	$(M*1000 *Z1) / (Z2*\delta)$
5.	MEQ / L	1000
6.	G / L	$(M * Z1) / Z2$
7.	MOLARITY(T)	$Z1 / (C*Z2)$
8.	—	1

2) Constants (12-20) for sample by weighing (mg.)

CODE NO.	UNITS	FORMULA
12.	MOLARITY(T)	$(M*Z1) / (1000 *Z2)$
13.	% ASSAY	$(M*Z1) / (Z2*10)$
14.	MG / G	$(M*Z1) / Z2$
15.	PPM	$(M*1000*Z1) / Z2$
16.	MOL / KG	$Z1 / Z2$
17.	ML / G	1
18.	TAN	56.1
19.	TBN	56.1
20.	—	1

- Z1 : Valency of the titrant. (Mol. wt. / Eq. Wt.)
Z2 : Valency of the sample. Eq. Wt. depends upon the nature of titration reaction.
M : Molar Mass (gms.).
C : Concentration of sample (mol / L).
TAN : Milligrams of KOH required to neutralize 1 gm. of sample.
TBN : Basicity equivalent to no. of mg. of KOH per gram of sample.
 δ : Density of sample

Note : Constant 8 is for getting no result unit printed for liquid sample and Constant 20 is for getting no result unit printed for solid sample.(Only volume is printed)

NOTE : Auto Incremented Run number is printed on the right side of every printout. If the volume reqd. is negative (blank/back volume exceeds titration volume) the result is printed and displayed **NEGATIVE**.

• Determination of Molarity:

Standardizing 0.1 M Sodium Hydroxide (NaOH) by Potassium Hydrogen Phthalate (KHP). Enter Constant related to Standard taken volumetrically in Code 7 or enter in Code 12 if Standard is taken by weight (Refer Page 71)

Constant Value = $(M \cdot Z1) / (1000 \cdot Z2) = (204.22 \cdot 1) / (1000 \cdot 1) = 0.2042$ in Code 12. Z1 is valency of the titrant and Z2 is valency of sample. Also enter the weight of sample taken for analysis.

- Select method for standardization.
- Perform titration with Reference Standard. (Refer Run Method Page).
- After titration Display reads the result.

Vol. Reqd : 0.9425 ml Result : 1.0025 Molarity (T)

Result – 1) Print 2) Reprocess 3) Display 4) Display Stat. 5) Run Last Method:

After repeat runs take statistics report, the mean value is displayed / printed. The mean Molarity (T) can be put in the required prog. no. manually. (for manual entry of molarity see below).

NOTE : Range for entering Molarity value is from 0.0001 to 99.99999 mol./lit.

MAINTENANCE AND CLEANING OF BURETTE/SYRINGE

- Periodic cleaning of the burette, plunger adapter and tubing with suitable solvent is recommended.
 - Do not dry the Non glass parts in oven.
 - Never place O-rings or rubber washers in organic solvents like chloroform, carbon tetra chloride.
1. Select the Burette option in Function Menu to remove the Assembly from the titrator unit.
 2. Follow the Display Screen instructions until the remove burette instruction is read on the display.
 3. Clean the syringe with suitable solvent. Be careful while handling plunger teflon tip against damaging as it may result in leakage.
 4. To mount it again, first locate the plunger of burette in the plunger drive arm and slide the assembly until click sound is heard with display message for type of assembly.
 5. Follow the Display Screen instructions until the system is restarted with correct burette assembly display message
 6. After Initialization and self-check is over perform rinsing of the entire liquid path – tubing, valve and burette.

To eliminate air bubbles from the burette, carry out at least 2-3 syringe strokes by PRIME.

MAINTENANCE AND CLEANING OF ELECTRODES

The glass electrode should be cleaned occasionally in moderately concentrated acid and kept in suitable buffer when not in use. See also the maintenance instructions packed with the electrode as set forth by the makers. The following procedure should be adopted for cleaning the glass electrode:

- 1) Rinse the electrode in pure water by changing water for more than 3-4 times, and wipe off water with clear filter paper or absorbent cotton.
- 2) If the glass electrode is extremely dirty, clean it with 0.1 N hydrochloric acid, soapsuds or detergent as necessary for a short time, and clean it thoroughly with pure water.
- 3) The electrode, which has been left dry for a long time, should be put in 4.00pH buffer to reach equilibrium before use.
- 4) If any foreign material is contained in porous frit (containing solution) of the reference electrode abnormal potential response may occur. Clean the solution area thoroughly with pure water.
- 5) After rinsing the electrode, clean their tips with pure water and wipe off water with clean filter paper or absorbent cotton. **DO NOT RUB THE SENSING PART.**
- 6) Now connect the electrode and observe the response of electrode as mentioned.

APPENDIX A


ALPHANUMERIC ENTRY

Alphanumeric entries are done to enter a Sample name, Titrant name (ENTER key) and Identification no. (RUN mode). Maximum 20 alphanumeric characters can be entered for any alphanumeric entry.

e.g. To enter Sample Name : say 2 HEXANE.

Please follow sequence of operations

NOTE: Press(2) means press the button twice.



Similarly enter Identification Number and Titrant Name.

APPENDIX B

ERROR CODES

# If the number of doses exceeds 200, during any titration run, titration run is aborted, error message is displayed and 'Run aborted' is printed on report.	ERROR 1
# If the titration run volume exceeds 99.99 ml, during any titration run, titration run is aborted, error message is displayed and 'Run aborted' is printed on report.	ERROR 1
# If after 2 doses Ctrl. Lt. for pH titration exceeds sample pH, error message is displayed and the titration run is aborted and 'Run aborted' is printed on report.	ERROR 1
# If the second element is 0 (i.e. this condition occurs only after resetting the system), graph cannot be printed, error message is displayed.	ERROR 2
# If number of doses exceeds 199, during any titration run, graph cannot be printed, error message is displayed.	ERROR 2
# If 1 ST deriv value is out of range (-32000 to 32000), during any run, graph cannot be printed, error message is displayed.	ERROR 3
# If the printer is not connected while printing data report or method parameter report, error message is displayed.	Printer not Ready
# If the method to be scanned or run is not already entered, error message is displayed.	Method does not exist

APPENDIX B

ERROR CODES (Contd.)

- | | |
|---|---|
| # If the electrode's offset or buffer values are not within specified range, error message is displayed. | Incorrect Buffer or Faulty Electrode |
| # If the method parameters are not suitable for running the method by changing burette, error message is displayed. | Check Method Parameter |
| # If the entered parameter is not within the specified range, error message is displayed. | Invalid Entry |
| # If the proper electrode is not connected, or BNC is open (i.e. electrode is not connected), error message is displayed. | Connect Proper Electrode |

APPENDIX C

PRINTER DIP SWITCH SETTINGS

The printer has twelve DIP switches grouped in two groups, mounted on the back panel. DIP switch 1-1 is the switch at the far left side, and the one at the far right is DIP switch 2-4.

The following table summarizes the settings of switches for WIPRO LX - 800 printer.

Switch Number	Position
SW 1 – 1	DOWN (OFF)
SW 1 – 2	UP (ON)
SW 1 – 3	UP (ON)
SW 1 – 4	DOWN (OFF)
SW 1 – 5	DOWN (OFF)
SW 1 – 6	UP (ON)
SW 1 – 7	UP (ON)
SW 1 – 8	UP (ON)
SW 2 – 1	UP (ON)
SW 2 – 2	DOWN (OFF)
SW 2 – 3	DOWN (OFF)
SW 2 – 4	DOWN (OFF)

Note : When a change is made in DIP switch setting, turn off the power, reset the switch or switches, then turn on the power again. The printer checks and recognizes new settings only at the time power is turned on.

**APPENDIX D
PARAMETER SELECTION TABLE**

PARAMETER	DESCRIPTION	DEFAULT	LIMITS
<i>Incremental Mode</i>			
PREDISPENSE VOLUME	Initial Predispense Volume	0	Maximum 95000 µl.
PREDISPENSE DOSE	Dosing size after each pre-dispense interval	OFF	Factor of PDV & In steps of 1/1000th of burette.
PREDISPENSE INTERVAL	Time between two doses during predispense	OFF	0 to 999 sec in step of 1 sec
STIRTIME	Time between predispense end and actual run start	'5'	0 to 999 sec in step of 1 sec
INITIAL DOSE	Dosing size after each interval	OFF	1/1000th of burette to 1/5th of burette in steps of 1/1000th of burette.
INTERVAL	Time between two doses	OFF	1 to 999 sec in step of 1 sec
VOLUME LIMIT	Maximum volume of titration	'10000'	> than initial dose and > than predispense volume Maximum 99.999 ml
END CRIT.	By VolLimit or EP	'EP'	9 EP or VL max.
<i>Equilibrium Mode</i>			
PREDISPENSE VOLUME	Initial Predispense Volume	0	Maximum 95000 µl.
PREDISPENSE DOSE	Dosing size after each pre-dispense interval	OFF	Factor of PDV & In steps of 1/1000th of burette.
PREDISPENSE INTERVAL	Time between two doses during predispense	OFF	0 to 999 sec in step of 1 sec
STIRTIME	Time between predispense end and actual run start	'5'	0 to 999 sec in step of 1 sec
INITIAL DOSE	Dosing size after each interval	OFF	1/1000th of burette to 1/5th of burette in steps of 1/1000th of burette.
EQ.FACTOR	Speed of titration	'5'	'1' to '99' in steps of 1
VOLUME LIMIT	Maximum volume of titration	'10000'	> than initial dose and > than predispense volume Maximum 99.999 ml.
END CRIT.	By VolLimit or EP	'EP'	9 EP or VL max.

APPENDIX D

PARAMETER SELECTION TABLE (Contd..)

PARAMETER	DESCRIPTION	DEFAULT	LIMITS
<i>Cutoff pH Mode</i>			
PREDISPENSE VOLUME	Initial Predispense Volume	0	Maximum 95000 µl.
PREDISPENSE DOSE	Dosing size after each pre-dispense interval	OFF	Factor of PDV & In steps of 1/1000th of burette.
PREDISPENSE INTERVAL	Time between two doses during predispense	OFF	0 to 999 sec in step of 1 sec
STIRTIME	Time between predispense end and actual run start	'5'	0 to 999 sec in step of 1 sec
INITIAL DOSE	Dosing size after each interval	OFF	1/1000th of burette to 1/5th of burette in steps of 1/1000th of burette.
INTERVAL	Time between two doses	OFF	1 to 999 sec in step of 1 sec
VOLUME LIMIT	Maximum volume of titration	'10000'	> than initial dose and > than predispense volume Maximum 99.999 ml.
END CRIT.	By VolLimit or EP	'EP'	9 EP or VL max.
CTRL. LIMIT	pH at which dose control starts	OFF	0-14 pH in steps of 0.01
END LIMIT	pH at which titration stops	OFF	0-14 pH in steps of 0.01
<i>pH stat</i>			
INITIAL DOSE	Dosing size after each interval	OFF	1/1000th of burette to 1/5th of burette in steps of 1/1000th of burette.
EQ.FACTOR	Speed of titration	'5'	'1' to '99' in steps of 1
VOLUME LIMIT	Maximum volume of titration	'10000'	> than initial dose and > than predispense volume Maximum 99.999 ml.
END CRIT.	By VolLimit or EP	'EP'	9 EP or VL max.
CTRL. LIMIT	pH at which dose control starts	OFF	0-14 pH in steps of 0.01
END LIMIT	pH at which titration stops	OFF	0-14 pH in steps of 0.01
RUN TIME	Maintaining pH within control band	OFF	1-9999 sec in steps of 1 sec.

APPENDIX E

TITRATION TYPE

- TYPE 1 : **ACID-BASE** In acid-base titration, an acid is determined by titration with a base, or vice-versa. The essential reaction is between H^+ and OH^- , giving water. The pH at the end-point depends on the dissociation constants of the reactants and products. It includes strong acid - strong base (e.g. HCl - NaOH), weak acid - weak base (e.g. KHP- Na_2CO_3), strong acid-weak base, weak acid-strong base and polybasic titrations.
- TYPE 2 : **NON-AQUEOUS** This includes the weak acids or bases which give more satisfactory end-points when titrations are carried out in non-aqueous media (organic solvents) than aqueous solutions. Because many of the solvents used are aggressive, volatile and obnoxious, non-aqueous titrations are carried out in closed environment. Compounds that may be determined by nonaqueous titrimetry include amines, amino acids, phenols, carbonyl compounds by oxidation and Schiff's bases.
- TYPE 3 : **REDOX** In these titrations a reducing agent is titrated with an oxidising agent and vice-versa. The common oxidising titrants are potassium permagnate ($KMnO_4$), potassium-di-chromate ($K_2Cr_2O_7$), iodine (I_2), potassium iodide (KIO_3), potassium bromate ($KBrO_3$) and sodium hypochlorite ($NaClO$). The most important reducing agents are ferrous ammonium sulphate ($(NH_4)_2Fe(SO_4)_2 \cdot 6H_2O$), sodium thiosulphate ($Na_2S_2O_3 \cdot H_2O$) and arsenic oxide(As_2O_3).
- TYPE 4 : **PRECIPITATION** These are titrations in which the analyte and titrant react to form a precipitate. The only common titrant used is silver nitrate (argentometric titrations) and its use is mainly restricted to the determination of chloride, bromide, iodide, cyanide and thiocyanate, which precipitates with silver ion.
- TYPE 5 : **COMPLEXIMETRIC/EDTA** Compleximetric titrations are used mainly to determine metal ions by use of complex forming reactions. In practise the titrants are compounds having iminodiacetic acid functional group, such as EDTA (ethylenediamine tetraacetic acid).
- TYPE 6 : **BACK TITRATION** : When the titration time is long and typical samples which do not react directly with titrant, i.e. need an intermediate are done by this method.
- TYPE 7 : **PHOTOMETRIC TITRATION** : Using Phototrodes this titration can be performed using optical principles.
- TYPE 8 : **THERMOMETRIC TITRATION** : Using temperature sensing principle this titration can be performed.

APPENDIX F

TITRATION MODES

- MODE 1: **INCREMENTAL** In this mode of titration fixed dose can be given at fixed interval of time up to the number of selected end point. This is used for determining sequential endpoints in multi-component sample. The titration is performable up to maximum 9 endpoints evaluation. The titration is stopped when volume limit is reached or number of number of EP set in method is reached.
- MODE 2: **EQUILIBRIUM** This is the most convenient and widely used mode. Dose as well as time gets adjusted depending on the trend of reaction and the rate of change of mV. The maximum dose during titration is the dose entered. Large mV change increases the time interval between the doses and the dose size reduces rapidly to minimum towards the end point. For slow reactions equilibrium factor should be raised. A smaller end dose is preferred for higher reproducibility and accuracy. This mode minimizes titration time while maintaining high precision. The data of titration curve are checked for the presence of an equivalence point. If an equivalence point is found and EP cutoff is selected the titration is stopped and the result is calculated. Only equivalence points with the first derivative greater than threshold are detected. Maximum 9 endpoints can be evaluated.
- MODE 3: **CUT-OFF BY pH** This is the same as equilibrium titration, but is based on pH values. In this mode selection of control and end limit is important, because it decides when the dose control begins during titration and up to what pH value the titration should be performed.
- MODE 4: **pH STAT** This is same as Cutoff by pH titration, since after attaining the set control pH, it starts giving controlled doses at regular intervals to maintain the pH for defined run time (Stat Run Time). The titration ends when set run time is reached. Volume consumed at particular time during the run can be calculated and used for result evaluation.

APPENDIX G

METHOD PARAMETERS

a) **INITIAL DOSE (μL)** :

This is dispense volume of titrant in μl . This volume gets added periodically during sample run. Dose selection should be such that significant change of mV value occurs after addition of every dose. Minimum dose that can be set is 5 μl . 1/5th of the burette volume is the maximum dose that can be entered in all modes of titration. The data storage can accommodate information for maximum 99 dosed increments; equivalence points can therefore only be recognized and evaluated if they lie within 99 increments.

b) **EQUILIBRIUM FACTOR** :

This factor changes the titration time for equilibrium and cutoff titration. The equilibrium titration time increases with increased factor value. Default value '2' is suitable for most of the titrations. However one can enter these values from 1-99 depending upon nature of reaction etc. for Equilibrium and Cut-off modes. For slow/weaker reactions increase the equilibrium factor value suitably.

c) **INTERVAL (sec)** :

This is the delay time between dose for Incremental end point titrations (1-999 sec). It is important for following reasons:

- 1) To take care of Response time of the electrode.
- 2) To ensure that dispensed reagent is sufficiently mixed before dispensing next dose.

d) **STAT INTERVAL (sec)** :

This is the time between recordings of the pH values during pH stat run. Start Interval for pH stat titration can be varied from 1 to 99 sec.

e) **CONTROL LIMIT (pH)** :

This is pH value at which the dose control occurs in Cut-off by pH titration. This should be chosen such that at least 5-6 doses are given before the end point detection. Control limit ranges from 0-14 pH for cutoff mode.

f) **END LIMIT (pH)** :

This is the pH value up to which cutoff by pH titration should be performed. If the end point is detected before the end limit, titration gets over not by the end limit but by the endpoint detection. End limit ranges from 0-14 pH for cutoff mode.

APPENDIX H

SAMPLE ANALYSIS PARAMETERS

a) **PRE-DISPENSE CRITERION (mV)** :

This is the mV value difference after which predispense will automatically stop and programmed doses will be given (0-1000). The value '0' allows the use of Pre dispense Volume function in which programmed doses are given at fixed minimum time of the system.

b) **PREDISPENSE VOLUME (μL)** :

In general, small predispensing is advisable since the electrode potential at the start of the titration changes drastically or is poorly defined. Predispense volume addition start the reaction without over-shooting equilibrium point and saving time. Pre-dispense should be such that at least 4-5 doses are given before endpoint. Max. predispense volume is 32 ml.

c) **PREDISPENSE DOSE (μL)** :

This is dispense volume of titrant in μl . This volume gets added periodically during predispense at the set predispense interval. Predispense dose should be multiple of 1/1000 th of burette and should add upto predispense volume. Burette volume is the maximum predispense dose that can be entered.

d) **PREDISPENSE INTERVAL (sec)**:

This is the delay time between predispense doses given in predispensing for all modes of titrations (0-999 sec).

e) **THRESHOLD** :

The threshold is principally used to prevent relatively small disturbances in the curve being mistaken for an equivalence point. The first derivative must be larger than threshold to ensure recognition of equivalence point. Default value '50' is suitable for most of the titrations. However one can enter these values from 1-32000 depending upon nature of titration for all modes. Additionally threshold can also be used to bypass smaller peaks.

f) **STIR TIME (sec)** :

Before measurement of the initial potential and start of the titration, a stirring time can be selected which allows homogenization and dissolution of the sample. Stir time ranges from 0-999 sec.

g) **VOLUME LIMIT (μL)** :

This is the volume up to, which the titration should be done, even if end point is reached and is applicable for all titration. Max. volume limit is 65 ml.

h) **MINIMUM CONTROL DOSE (μL)** :

This is the minimum dose that can be given at endpoint for equilibrium and cutoff by pH titration mode. Values enterable from 5-500 μL

APPENDIX H

SAMPLE ANALYSIS PARAMETERS (contd.)

- i) **STAT RUN TIME (sec)** :
This is the time for which pH-Stat titration should be continued after attaining set pH value. Run time can be entered for the range of 1-9999 sec.
- j) **VTIME 1 (sec)** :
This is the time at which one wants to calculate the volume consumed during pH stat titrations. VTime 1 can be entered in the range of 1-9999 sec. The value of VTime 1 should be multiple of stat interval time.
- k) **VTIME 2 (sec)** :
This is the time at which one wants to calculate the volume consumed during pH stat titrations. Vtime 2 can be entered in the range of 1-9999 sec. The value should be higher than VTime 1 and should be multiple of stat interval time.
- l) **FACTOR** :
This is a factor used for calculation of result (enzymatic activity) for pH stat titration. The value can be entered in the range 0.0001 to 10000.
- m) **BUFFER VALUE (mV)** :
This is the buffer value entered for Method 38 and 40 i.e. for ASTM D664 & D4739 respectively. The endpoint can be calculated for the entered buffer value, the titration ends at 100 mV above the set buffer value. Value can be entered in the range 1-3200 mV.
- n) **BUFFER SIGN** :
This is the polarity for the buffer value entered for Method 38 and 40 i.e. for ASTM D664 & D4739 respectively as mentioned above. The endpoint can be calculated for the entered buffer value, the titration ends at 100 mV above the set buffer value. + or - sign can be assigned to the entered buffer value.
- o) **TITRATION TREND** :
This is the decided by the titration data table. If mV rises with every dose given i.e. (more positive with every dosing or negative to positive with every dosing) the trend is + (positive). If mV falls or drops with every dose given i.e. (more negative with every dosing or positive to negative with every dosing) the trend is - (negative). This is again related to ASTM D664 & D4739 i.e. method 38 & 40 run.

APPENDIX I

CALCULATION PARAMETERS

- a) **CONCENTRATION** :
This is the molarity of titrant to be considered for result calculation. For standardization of titrant, before run the value entered is 1 and then evaluated value gets transferred in respective code. Concentration value ranges from 0.0001-99.9999 mol/lit.
- b) **PRIMARY DOSE** :
This is volume in ml, which gets added in final volume required for titration. Prime the required volume, enter the prime volume as primary dose before running the required program. After running the program primary dose entered will get added to total volume. Primary dose ranges from 0.0001-99.9999 ml.
- c) **BLANK VOLUME** :
This is the volume in ml, which gets subtracted in final volume required for titration. This is determination of titrant consumed for neutralization by equivalent volume of solvent to be used along with the sample. Blank volume ranges from 0.0001-99.9999 ml.
- d) **CONSTANT** :
This is the constant related to the sample to be entered for result calculation. Constant code 1-10 are for liquid samples and 11-20 are for solid samples. Constant value ranges from 0.0001-999999.99.
- e) **SAMPLE VOLUME** :
The amount of sample taken for run to be considered for result calculation for liquid samples. Sample volume ranges from 0.0001-99.9999 ml.
- f) **SAMPLE WEIGHT** :
The amount of sample taken for run to be considered for result calculation for solid samples. Sample weight ranges from 0.0001-99.9999 gms.
- g) **VOLUME FACTOR** :
This is factor which gets multiplied to the blank volume entered for back titration. It is used for calculation purpose. Volume factor can be entered in the range 0.1000-5.0000.

APPENDIX J

Standard Methods (41-50) :

Method 41 : Standardizing 0.1 M NaOH against KHP

METHOD NO	: 41	TITRANT CONC	: 0.1 M (appx.)
SAMPLE NAME	: KHP	TITRANT NAME	: NAOH
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 100 mg (apx)
TITRATION TYPE	: Acid-Base Aq.	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 3600 µl
PRE.DISP.DOSE	: 200 µl	PRE.DISP.INTV.	: 2 sec
STIR TIME	: 30 sec	INITIAL DOSE	: 100 µl
EQ. FACTOR	: 6		
VOLUME LIMIT	: 8000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 300	MIN. CONTROL DOSE	: 10 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.2042	CONSTANT CODE	: 12
REPORT	: None		

Method 42 : Standardizing 0.1 M HCl against

METHOD NO	: 42	TITRANT CONC	: 0.1 M (appx.)
SAMPLE NAME	: Na ₂ CO ₃	TITRANT NAME	: HCL
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 35-40mg (ap)
TITRATION TYPE	: Acid-Base Aq.	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 5000 µl
PRE.DISP.DOSE	: 200 µl	PRE.DISP.INTV.	: 1 sec
STIR TIME	: 30 sec		
INITIAL DOSE	: 100 µl	EQ. FACTOR	: 6
VOLUME LIMIT	: 10000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 100	MIN. CONTROL DOSE	: 5 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.0530	CONSTANT CODE	: 12
REPORT	: None		

APPENDIX J
Standard Methods (41-50) contd ... :

Method 43 : Standardizing 0.1 M Alc. KOH against KHP

METHOD NO	: 43	TITRANT CONC	: 0.1 M (appx.)
SAMPLE NAME	: KHP	TITRANT NAME	: Alc. KOH
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 100 mg (apx)
TITRATION TYPE	: Acid-Base Aq.	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 3600 µl
PRE.DISP.DOSE	: 200 µl	PRE.DISP.INTV	: 1 sec
STIR TIME	: 30 sec		
INITIAL DOSE	: 100 µl	EQ. FACTOR	: 6
VOLUME LIMIT	: 8000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 100	MIN. CONTROL DOSE	: 10 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.2042	CONSTANT CODE	: 12
REPORT	: None		

Method 44 : Standardizing 0.1 M HClO₄ against KHP

METHOD NO	: 44	TITRANT CONC	: 0.1 M (appx.)
SAMPLE NAME	: KHP	TITRANT NAME	: HClO ₄
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 100 mg (app)
TITRATION TYPE	: Non Aqueous	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 3600 µl
PRE.DISP.DOSE	: 200 µl	PRE.DISP.INTV.	: 2 sec
STIR TIME	: 30 sec		
INITIAL DOSE	: 100 µl	EQ. FACTOR	: 6
VOLUME LIMIT	: 8000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 300	MIN. CONTROL DOSE	: 10 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.2042	CONSTANT CODE	: 12
REPORT	: None		

APPENDIX J

Standard Methods (41-50) contd :

Method 45 : Standardizing 0.1 M NaNO₂ against Sulphanilic Acid

METHOD NO	: 45	TITRANT CONC	: 0.1 M (appx.)
SAMPLE NAME	: SULPHANILIC ACID	TITRANT NAME	: NANO ₂
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 100 mg (apx)
TITRATION TYPE	: Redox	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 4000 µl
PRE.DISP.DOSE	: 200 µl	PRE.DISP.INTV.	: 5 sec
STIR TIME	: 30 sec		
INITIAL DOSE	: 100 µl	EQ. FACTOR	: 10
VOLUME LIMIT	: 7000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 300	MIN. CONTROL DOSE	: 10 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.1732	CONSTANT CODE	: 12
REPORT	: None		

Method 46 : Standardizing 0.02 M KMnO₄ against Na Oxalate

METHOD NO	: 46	TITRANT CONC	: 0.02M(appx.)
SAMPLE NAME	: NA OXALATE	TITRANT NAME	: KMNO ₄
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 45-50mg (ap)
TITRATION TYPE	: Redox	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 6000 µl
PRE.DISP.DOSE	: 100 µl	PRE.DISP.INTV.	: 5 sec
STIR TIME	: 10 sec		
INITIAL DOSE	: 100 µl	EQ. FACTOR	: 15
VOLUME LIMIT	: 10000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 500	MIN. CONTROL DOSE	: 10 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.3350	CONSTANT CODE	: 12
REPORT	: None		

APPENDIX J

Standard Methods (41-50) contd ... :

Method 47 : Standardizing 0.1 M Na₂S₂O₃ against KIO₃

METHOD NO	: 47	TITRANT CONC	: 0.1 M (appx.)
SAMPLE NAME	: KIO ₃	TITRANT NAME	: NA ₂ S ₂ O ₃
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 25-30 mg (ap)
TITRATION TYPE	: Redox	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 7000 µl
PRE.DISP.DOSE	: 200 µl	PRE.DISP.INTV.	: 1 sec
STIR TIME	: 30 sec		
INITIAL DOSE	: 100 µl	EQ. FACTOR	: 5
VOLUME LIMIT	: 10000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 300	MIN. CONTROL DOSE	: 5 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.0356	CONSTANT CODE	: 12
REPORT	: None		

Method 48 : Standardizing 0.02 M Hg(NO₃)₂ against NaCl

METHOD NO	: 48	TITRANT CONC	: 0.02M(appx.)
SAMPLE NAME	: NaCl	TITRANT NAME	: Hg(NO ₃) ₂
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 15mg (ap)
TITRATION TYPE	: Redox	TITRATION MODE	: Incremental
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 5000 µl
PRE.DISP.DOSE	: 500 µl	PRE.DISP.INTV.	: 2 sec
STIR TIME	: 30 sec		
INITIAL DOSE	: 50 µl	INTERVAL	: 5 sec
VOLUME LIMIT	: 8000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 10		
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.1167	CONSTANT CODE	: 12
REPORT	: None		

APPENDIX J

Standard Methods (41-50) contd :

Method 49 : Standardizing 0.1 M AgNO₃ against NaCl

METHOD NO	: 49	TITRANT CONC	: 0.1 M (appx.)
SAMPLE NAME	: NaCl	TITRANT NAME	: AgNO ₃
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 40 mg (apx)
TITRATION TYPE	: Precipitation	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 5000 µl
PRE.DISP.DOSE	: 100 µl	PRE.DISP.INTV.	: 1 sec
STIR TIME	: 30 sec		
INITIAL DOSE	: 100 µl	EQ. FACTOR	: 6
VOLUME LIMIT	: 10000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 100	MIN. CONTROL DOSE	: 10 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.0584	CONSTANT CODE	: 12
REPORT	: None		

Method 50 : Standardizing 0.1 M EDTA against ZnSO₄

METHOD NO	: 50	TITRANT CONC	: 0.1 M (appx.)
SAMPLE NAME	: ZnSO ₄	TITRANT NAME	: EDTA
MEASURE. TYPE	: Potentiometric	SAMPLE WEIGHT	: 100 mg (app)
TITRATION TYPE	: Comp./EDTA	TITRATION MODE	: Equilibrium
PRE.DISP.CRIT	: 0 mV	PRE.DISP.VOLUME	: 0 µl
PRE.DISP.DOSE	: 0 µl	PRE.DISP.INTV.	: 0 sec
STIR TIME	: 0 sec		
INITIAL DOSE	: 200 µl	EQ. FACTOR	: 5
VOLUME LIMIT	: 5000 µl	END CRITERIA	: EndPoint (1)
THRESHOLD	: 10	MIN. CONTROL DOSE	: 10 µl
MOLARITY (mol/lit.)	: 1.0000		
BV/PD/NONE	: None		
CONSTANT	: 0.2875	CONSTANT CODE	: 12
REPORT	: None		

APPENDIX K

K905XX – Data Downloading Capability to Personal Computer

Installation/setup Instructions:


The instrument software has been incorporated with RS232C protocol interface for transferring the printable data to PC having standard data terminal software Windows - **Hyper terminal**.

Features of PC data transfer,

1. Data transfer selectivity in function menu : Printer or Personal computer.
2. All reports - a) Result Report b) Data c) Condensed Report d) Statistical Data e) Calibration Report of pH f) Methods which are available for printing will be transfer to PC (needs Hyper terminal active in the PC).
3. The data is captured as ASCII character via Hyper terminal software and then can be saved by the user as MS word doc or text file by Note Pad with suitable name by using Cut Paste utility.

COM Port Configuration/Settings:-


Setting HYPERACCESS NEW CONNECTION and COM port selection:-






If the COM port chosen is not free or being used by another application a pop up screen will be displayed with message “Unable to Open COM”.

When a free Com port is selected The COM port Properties Pop up will be displayed on the screen.






Set the Comm Port Properties, Baud-9600, Data bits 8, Parity none, Stop bit 1, Flow Control- None or Xon/Xoff.



➤ Once the setting are done Save the Session so as the Session file is generated and can be run again to download the data from **PC-TITRA**.

Now From the terminal window FILE menu select the Properties tool to set the Terminal Emulation mode and Data Line
 Once You click on the Properties menu - Properties window will appear.



Select Settings option to set the Emulation mode and function setting.


Function:- Arrow & Ctrl Key as Terminal key.

Backspace key sends: Ctrl+H

Emulation :- VT100


Telnet id:- VT100.

Terminal Setup:- Character Set = ASCII.



Now Setup the ASCII parameters

Click the ASCII setup icon to enter the setting window.



- Click the Boxes as shown in the graphics for**
- Line end with line feed,
 - Echo typed characters,
 - Append line feed to incoming line ends,
 - wrap lines that exceed terminal width,
 - Then Click OK to return to Properties window

Click OK Icon on the Properties window to complete the settings.

Now once again Save the Session to store the settings.

To Generate the Icon for the TITRA connectivity with PC, right click on the file name being in the Hyper terminal Window – select Create Shortcut to desk top.

Click on the Icon and start the application. An Terminal Window will open automatically and now the PC is ready to receive the data from K905XX - Automatic Potentiometric Titrator.