



designed for scientists



IPCS Pipette Calibration Software

/// Data Sheet

To ensure the accuracy and precision of each pipette, regular calibration checks are necessary. IPCS provides user-friendly pipette calibration by scheduling, performing and documenting all calibration checks. Perform and schedule your next calibration with ease. Automatic transfer of measurements from the weighing balance helps improving the quality of calibration. The software comes equipped with a special auto mode which helps in detecting and transmitting temperature, humidity and pressure.

- Calibrate your pipette on your own
- Features all variations of IKA Pipettes

www.ika.com

Subject to technical changes





designed for scientists

- User-friendly interface
- Easy scheduling of next calibration
- Automatic transfer of readings from the weighing balance
- Ease of selecting standards for calibration checks
- Store your comments and results in the calibration log
- Automatic calculation of the Z factor





designed for scientists

