



KPM-203D-2 / Polishing Machine / Preparing Specimens

### Summary

The tester is desgned to obtain specular surfaces of various specimen It is useful for observing the specular surfaces obtained by a grinder w microscope.

#### Features

This machine is desgned or manufactured to obtain specular surfaces specimens, such as metalic materials.

#### DETAILED INFORMATION

## **▶** Specification

■ Number of rotation: 5 ~ 900rpm knob controller

Disc Size : AL Φ203 ~ 254 mm, tapered type

Cooling Water: Free-form valve mounted

Cleaning method by water spray in bowl

Motor system for minimizing noise and vibration

SOL mountable for cooling water supply and shut off

• Economical price and new concept of consumables included

• 3 automatic head : Optional

• Size : 650(W) \* 500(D) \* 400(H) mm

# **▶** Application

Molded specimen, products that require gloss