

# SL 3111 Validation-Set for validation of large steam sterilizers

## for large steam sterilizers

For the validation of steam sterilizers according to ISO 17665 as well as for the validation of washer-disinfectors and washer-disinfectors for endoscopes according to ISO 15883. The set can individually be expanded or compiled yourself from one or more data loggers (EBI 12 or EBI11), the appropriate interface and the software Winlog.validation.  
Film: [download](#)



## Product Description

The set contains:

- 5 x EBI 12-T441 Temperature data loggers
- 1 x EBI 12-TP226 Temperature/ pressure data logger with AL 101 silicone protection box
- EBI IF 200, 4-port Interface with USB connection and antenna
- Software Winlog.validation
- Carrying case „SYSTAINER“

Part No.: 1250-3101

## Technical data

-----

---

## Software

### Winlog.pro Software

Winlog.pro - The general software for programming and evaluation of ebro data loggers

Part No. 1340-2355



Hardware requirements	
	Processor: Pentium 1 GHz, Main memory: 1 GB, Free hard disk space: 1 GB, USB

<b>Software requirements</b>	Windows Vista (32 bit or 64 bit), Windows 7,8,10 (32 bit or 64 bit)
<b>Further requirements</b>	Resolution: 1024x768

---

## Winlog.med Software

Software (Standard-Software for the CSSD)



Part No. 1340-2363

<b>Hardware requirements</b>	Processor: Pentium 1 GHz, Main memory: 1 GB, Free hard disk space: 1 GB, USB
<b>Software requirements</b>	Windows Vista (32 bit or 64 bit), Windows 7,8,10 (32 bit or 64 bit)

---