


HI931



Thoughtfully designed. Inside and out.




This new generation of titrator features a 50% smaller footprint for maximum use of your lab space. The sealed body is precision-engineered to resist spills and dust and a high-contrast display shows you everything you need to know during testing.

With an angled electrode holder and removable stirrer, smaller volume titrations are now seamless. Not only can your titrator hold more electrodes and tubes, but it is also easier than ever to reposition your electrode holder using a newly designed press-to-release button and rotating disk.

Our unmatched 40,000-step dosing pump is capable of dosing extremely small volumes of titrant to help you achieve a very precise endpoint for greater consistency. There is no other like it. Add a second analog board, burette, or peristaltic pump for more power, more usability, and more versatility.

Our exchangeable Clip Lock™ burette system helps you conveniently separate and switch between your reagents so you don't have to worry about cross-contamination.

A dynamic dosing algorithm allows for timely and accurate titration results by providing larger doses near the beginning of a titration and smaller, more precise doses near the titration endpoint.



Set a new standard for titration.

The HI931 titrator delivers incredible results and speed for increased productivity. Titrate for a variety of measurements including acids, bases, redox, and selective ions. With no additional programming upgrades to purchase, you can start measuring right away.

Ordering
Information

HI931-01 (US plug (type A)) and **HI931-02**: (European plug (type C)) titrator with one analog board, propeller stirrer with stand, 25 mL glass burette, dosing pump, temperature sensor, USB cable, USB flash drive and PC software.



www.hanna.co.za

JHB - 011 615 6076
DBN - 031 701 2711
CT - 021 946 1722

HANNA
instruments

HI932



Thoughtfully designed. Inside and out.




This new generation of titrator features a 50% smaller footprint for maximum use of your lab space. The sealed body is precision-engineered to resist spills and dust and a high-contrast display shows you everything you need to know during testing.

With an angled electrode holder and removable stirrer, smaller volume titrations are now seamless. Not only can your titrator hold more electrodes and tubes, but it is also easier than ever to reposition your electrode holder using a newly designed press-to-release button and rotating disk.

Our unmatched 40,000-step dosing pump is capable of dosing extremely small volumes of titrant to help you achieve a very precise endpoint for greater consistency. There is no other like it. With the second analog board, a burette or peristaltic pump can be added for more power, more usability, and more versatility.

Our exchangeable Clip Lock™ burette system helps you conveniently separate and switch between your reagents so you don't have to worry about cross-contamination.

A dynamic dosing algorithm allows for timely and accurate titration results by providing larger doses near the beginning of a titration and smaller, more precise doses near the titration endpoint.



Set a new standard for titration.

The HI932 titrator delivers incredible results and speed for increased productivity. Titrate for a variety of measurements including acids, bases, redox, and selective ions. With no additional programming upgrades to purchase, you can start measuring right away.

Ordering
Information

HI932-01 (US plug (type A)) and **HI932-02**: (European plug (type C)) titrator with two analog boards, propeller stirrer with stand, 25 mL glass burette, dosing pump, temperature sensor, USB cable, USB flash drive and PC software.



www.hanna.co.za

JHB - 011 615 6076
CT - 021 946 1722
DBN - 031 701 2711

HANNA[®]
instruments