

TRANSCUTANEOUS ELECTRICAL NERVE STIMULATION (TENS) TRAINER**Model Number : GOTT-TENS-01****DESCRIPTION**

GOTT-TENS-01 is designed to demonstrate the operation, theory, repair and servicing of Transcutaneous Electrical Nerve Stimulation (TENS) equipment that used at clinics & hospitals. The unit is given the fault diagnosis so as to facilitate the teaching of trouble shooting techniques. A fault simulator is provided with circuit block diagram of the device for the students to study theoretically and experimentally in order to understand the operation of the device.

The fault is activated by button that disables or create malfunction of certain circuit or components so that the device become faulty.

Based on the symptoms and the signal at the test points, the students then do the troubleshooting exercise to find the faulty circuit or components. After locating the failure, the fault can be reset, and then the device work normally.

FEATURES

- Real Transcutaneous Electrical Nerve Stimulation unit
- High quality industrially accepted component
- Bold color identification of circuits
- Designed to simulate common faults
- Detailed instruction & experiment manual

SPECIFICATIONS

- 4.3 inch color LCD with backlight
- Three working modes: Rate/Volume/Time mode
- Heating and remote control functions
- Audible and Visible alarms
- Circuit Block Diagram of Transcutaneous Electrical Nerve Stimulation
- Electrical Test Points
- Fault Simulation

EXPERIMENT TOPICS

- Understanding Transcutaneous Electrical Nerve Stimulation Unit
- Block Diagram Transcutaneous Electrical Nerve Stimulation Explanation
- Operating The Transcutaneous Electrical Nerve Stimulation Training Unit
- Troubleshooting of Rectifier
- Troubleshooting of Microcontroller
- Troubleshooting of the Battery
- Troubleshooting of Control Panel
- Troubleshooting of Signal Conditioning

Manuals :

- (1) All manuals are written in English
- (2) Model Answer
- (3) Teaching Manuals

General Terms :

- (1) Accessories will be provided where applicable.
- (2) Manual & Training will be provided where applicable.
- (3) Design & specifications are subject to change without notice.
- (4) We reserve the right to discontinue the manufacturing of any product.

Warranty :

2 years

ORDERING INFORMATION :

ITEM	MODEL NUMBER	CODE
TRANSCUTANEOUS ELECTRICAL NERVE STIMULATION (TENS) TRAINER	GOTT-TENS-01	200-409

* Proposed design only, subject to changes without any notice