

“New-P.NIX® UP 900” Ultra Pure(UP) Water Purification System, Max 1.5L/min
 With Total Organic Carbon(TOC), 2-Steps of Filter Exchange Indicator, (UP) Up to 18.3MΩ·cm, 초순수 제조장치



- Certified CE, ISO 9001, ISO 14001
- RS 232C Port and Volumetric Dispensing Function by Timer
- Auto Regular Recycling to Maintain the Optimum water Quality Level
- Programming by 2 Steps for the Filter Exchange
- Filter Replace Message Indicating by LCD Display & Blinking LED Light

- * Microprocess One Touch Technology Controller
- * 16 × 4 Characters Backlight LCD Digital Display
- * Obtain GD, QD Mark, Awarded the Korean Excellent Products
- * Included 0.2μm Final Capsule Filter, UP Membrane Pack
- * Prerequisite : Pure Water is Required (DH.WWP.RO180)

Specifications & Ordering information

Spec.	Model	New-P.NIX UP 900			
		Scholar-type	UV-type	Bio-type	Integrate-type
Certified Mark		CE, ISO 9001, ISO 14001			
Product Flow Rate(UP)		Max. 1.5L/min			
Water Quality Display		• UP Product : Up to 18.3MΩ·cm • Feed water : 0 to 250μS/cm • Temperature : 0 to 99.9℃			
Display		16 × 4 Characters Backlight LCD Digital Display			
Included Filter Composition		• 0.2μm Final Capsule Filter • UP Pack (DH.WWP.UPS)	• 0.2μm Final Capsule Filter • UP Pack (DH.WWP.UPU) • 254/185UV Lamp	• 0.2μm Final Capsule Filter • UP Pack (DH.WWP.UPS) • MW 5000 UF Filter(L)	• 0.2μm Final Capsule Filter • UP Pack(DH.WWP.UPI) • 254/185UV Lamp • MW 5000 UF Filter(L)
Controller		Microprocess One Touch Technology			
Dimensions (w × d × h)		300 × 500 × h 470mm			
Power		210~230V, 50/60 Hz			
Cat. No.		DH.WWP.UP.S	DH.WWP.UP.U	DH.WWP.UP.B	DH.WWP.UP.I



(1) White, HDPE Aspiration Bottle

(2) Amber, HDPE Aspiration Bottle



(3) 0.2μm Final Capsule Filter



(4) P.NIX UP Pack, Scholar / Bio



(5) P.NIX UP Pack, Scholar-UV / Integrate

- TOC monitoring system is also available for the Integrate-type on request. (TOC monitoring display range : 0~200ppb)

Optional Accessory

SL.Bot3102	(1) Bottle, HDPE, Aspiration, White, 5L, Neck id-φ 70, φ 200 × h360mm
SL.Bot3103	, 10L, -φ 70, φ 250 × h460mm
SL.Bot3104	, 20L, -φ 95, φ 300 × h555mm
SL.Bot3105	, 50L, -φ 95, φ 380 × h700mm
SL.Bot3106	(2) Bottle, HDPE, Aspiration, Amber, 5L, Neck id-φ 70, φ 200 × h360mm
SL.Bot3107	, 10L, -φ 70, φ 250 × h460mm
SL.Bot3108	, 20L, -φ 95, φ 300 × h555mm
SL.Bot3109	, 50L, -φ 95, φ 380 × h700mm

※ Aspiration Bottles Details for Order : Please See “Bottles, Aspiration” Section.

Spare Filter & UV Lamp

Cat. No.	Description	Filter Replace Period
DH.WWP.0.2CF	(3) 0.2μm Final Capsule Filter	Approx. 6~8 Month
DH.WWP.UPS	(4) P.NIX UP Pack, Scholar / Bio	Approx. 1 Year
DH.WWP.UPU	(5) P.NIX UP Pack, Scholar-UV / Integrate	
DH.WWP.MWL	MW 5000 UF Filter(L)	Approx. 6~8 Month
DH.WWP.UV	254/185 UV Lamp	Approx. 3 Year

- The above mentioned Filter replacement period can be changed according to various factors such as feed water condition, producing volume and so on.