



Home (/) >> PRODUCT (<http://www.hssdgroup.com/PRODUCT/>) >> Metallographic Equipment (http://www.hssdgroup.com/PRODUCT/Metallographic_Equipment/)



P-2 metallographic sample polishing machine (combined type)

Product description:

1. Application: Model P-2 Metallographic specimen polishing machine is suitable for the polishing of the specimen which has been grinded. The specimen surface is very smooth after processing and can b



([http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/](http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/Compression_Tester/PC-0100/)) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/) (http://www.hssdgroup.com/PRODUCT/Compression_Tester/PC-0100/)

Share With:

We're here to help:
Easy ways to get the answers you need.

Send E-mail (<mailto:Louis@hssdtest.com>)

Chat Now (<skype:sales02@hssdtest.com?chat>)



(/)

Model P-2 Metallographic specimen polishing machine is suitable for the polishing of the specimen which has been grinded. The specimen surface is very smooth after processing and can be used to observe and measure the metallographic structure of specimen under microscope.

2. Main specifications:

Model	P-2
Polishing disc diameter	203mm
Rotating speed	1400r/min (special order)
power supply	AC 220v,50Hz (special order)
Net weight	30.75kg
Motor	180w
Dimension	888 × 478 × 915 mm