

Provide A Full Set Of Solutions For Battery Machines.



## 0.5mL To 50mL Laboratory Pipette StepMate Stepper

0.5mL to 50mL Laboratory StepMate Stepper

- Brand:TOB NEW ENERGY
- Item No.:TOB-StepMate
- Order(Moq):1set
- Payment:L/C,T/T,Western Union,
   Paypal
- Product Origin:China
- Shipping Port:XIAMEN
- Lead Time:20days

## **Product Detail**

0.5mL to 50mL Laboratory Pipette StepMate Stepper with seven sizes of disposable

## **SPECIFICATIONS**

- Light and ergonomic solution provides easy usability
- Single-handed operation
- Volume range from 10µL to 5mL
- Up to 48 dispensing steps
- Maintenance free
- Equipped with durable tip
- insertion lever
- Fits seven sizes of disposable polypropylene syringes from 0.5mL to 50mL

Dial setting	1	2	3	4	5
Number of steps	48	23	15	11	8

Cat. No.	Syringe volume (mL)	Sample volume (μL)				L)	A %	P %	Description
1790003	0.50	10	20	30	40	50	±0. 8	≤ 0.7	0.5mL, piston PE-HD, cylinder PP,1 packing,
									5 pcs/packing
1790003	1.25	25	50	75	100	125	± 0.8	≤0. 5	1.25mL, piston PE-HD, cylinder PP,1 packing,
									5 pcs/packing
1790003 4	2.50	50	100	150	200	250	± 0.8	≤0. 4	2.5mL, piston PE-HD, cylinder PP, 1 packing,
									5 pcs/packing
1790003 5	5.00	100	200	300	400	500	±0.	≤ 0.3	5mL, piston PE-HD, cylinder PP, 1 packing,
									5 pcs/packing
1790003	12.50	250	500	750	100	125 0	± 0.3	≤ 0.3	12.5ml, piston PE-HD, cylinder PP, 1 packing,
									5 pcs/packing
1790003	25.00	500	100	150	200	250	± 0.2	≤ 0.3	25mL, piston PE-HD, cylinder PP, 1 packing, 2pcs/packing with 1 pc adapter

1790003	50.00	100	200	300	400	500	±	≤	50mL, piston PE-HD,
8		0	0	0	0	0	0.2	0.2	cylinder PP, 1 packing,
									2pcs/packing with
									4 no adentes
									1 pc adapter

## PRODUCT DISPLAY

