



Bluepard Instruments Co., Ltd.
Suzhou Being Medical Device Co., Ltd.

Plant Growth Chamber

Microprocessor controller (with timing function)

Features:

- Microprocessor controller for temperature and humidity and light intensity, 30 programs
- Large LCD screen
- Simulate changing temperature and light in day/night
- Polished stainless-steel chamber, semicircular arcs at corners for easy cleaning, and the space between the shelves in the chamber is adjustable.
- R134a refrigerant, two imported compressors, continuous running time
- Independent temperature-limiting alarm system ensures experiments run safely.(option)
- RS485 connector is option which can connect computer to record the parameters and the variations of temperature.(option)



Specifications

Parameter	Model	MGC-250BP-2	MGC-350BP-2 MGC-450BP-2	MGC-350HP-2 MGC-450HP-2	MGC-800BP-2 MGC-1000BP-2	MGC-800HP-2 MGC-1000HP-2
Chamber Volume		250L	300L/450L		800L/1000L	
Temperature Range		With Lighting:10~50°C Without lighting:0~50°C				
Display Resolution		0.1°C				
Temperature Stability		±1°C				
Temperature Uniformity		±1.5°C			±2°C	
Humidity Range		BP-2: no humiditycontroller HP-2: 50~90%RH				
Humidity Accuracy		BP-2: no humiditycontroller HP-2: ±5~7%RH				
Lighting Intensity		0~12000LX	0~20000LX/0~25000LX		0~30000LX/0~35000LX	
Lighting Type		Shelves illumination	Three surface illumination		Shelves illumination(two Shelves)	
Programmer Function		Temperature, humidity, light intensity, can be set separately, and 30 programmers Time range is 1-59 hours for each programmer				
Power Consumption		1700W	1700W/2000W	1700W/2050W	3700W/3800W	
Electrical Requirement		220V 50HZ			380V 50HZ	
Ambient Temperature		RT+5~30°C				
Interior Dimension(W×D×H) mm		580×510×835	520×550×1140/700×550×1140		965×580×1430/900×580×1600	
Exterior Dimension(W×D×H)mm		725×740×1550	830×850×1850/950×850×1850		1475×890×1780/1410×890×1950	
Shelves		3(PCS)				

Options:

- CO₂ inlet
- CO₂ controller (imported IR CO₂ sensor)

•RS485 connector

