

## Baby Carriage Brake Performance Testing Machine TC-009



### Detailed Description

Made by CNS Hand-push Baby Carriage's effects standard, testing method is putting the baby model in the car, placed the car on the table with 60#sandpaper, look at the car if slipped when adjust the tilt angle of testing table 30 seconds later.

### Specifications:

Testing table's tilt angle: 0°-25° Adjustable

Surface: Put on 60# sandpaper

Area: 120cmx160cm

Dimension: 130x170x100cm

Adjust method: By hand