



Brush Pilling Tester GT-C20



Application:

ASTM D 3511 **Brush Pilling Tester** GT-C20 is used to determine the fabric pilling performance in simulate wearing and normal operating conditions, can be used for knitted and woven fabrics which used in cloth, home textile and automotive interiors decoration materials. It needs to be grade by a rating endoscopic , light cabinet or similar, and piling photographs.

Standards:

FORD BN-108-03, GM 9652P, ASTM D 3511

Key Specifications

Model	GT-C20
Circular sample holder	660g, soft foam rubber surface
Pilling range	19mm
Height of nylon brush	25.4mm
Square sample	Holder x6 pcs, rubber edge clamp, sandpaper backface
Rotate speed	58 r/min
Sponge size	51 x 102 x 152 mm
Power supply	220 V/50 Hz
Dimensions	914 x 305 x 381mm (L x W x H)
Weight	50kg
Standards	FORD BN-108-03, GM 9652P, ASTM D 3511
Standard accessories	Power line -1pc
Optional accessories	GT-C48 <u>Pneumatic Sample Press</u>